

# MAAT

medac Asparaginase Activity Test

For a reliable Control of Asparagine Depletion by specific Drug Monitoring



medac


## The medac Asparaginase Activity Test offers an Opportunity to optimize Asparaginase Treatment by therapeutic Drug Monitoring

- Monitoring of asparaginase activity can detect silent inactivation caused by the patient's immune response.
- Use of an alternative preparation in case silent inactivation is detected.
- Controlling the duration of asparaginase activity during treatment course.
- A target level of  $\geq 100$  U/L asparaginase activity is reported to ensure therapeutic benefit <sup>(4, 6, 7)</sup>.
- Determination of the lowest possible dosage of the drug, while retaining full efficacy.
- Rapid availability of exact results for an appropriate therapeutic decision.

### Kit Content

- Microplate, 6 x 8 wells (standard 96 well format), breakable, U-form
- 2 vials asparaginase control, app. 100 U/L, 0.5 ml, lyophilized
- 2 vials asparaginase calibrator, 0 U/L, 0.5 ml, lyophilized
- 2 vials asparaginase calibrator, 50 U/L, 0.5 ml, lyophilized
- 2 vials asparaginase calibrator, 300 U/L, 0.5 ml, lyophilized
- 2 vials asparaginase calibrator, 600 U/L, 0.5 ml, lyophilized
- 1 vial sample diluent, 110 ml, ready to use
- 2 vials substrate, 0.75 ml, ready to use
- 2 vials chromogen, 1 ml, concentrate
- 2 vials chromogen diluent, 2 ml, ready-to-use

## Test Principle




## Test Characteristics


- Stability: 8 weeks after opening the kit
- Analytic range: 30-600 U/L ( $\infty$  by further dilution, linearity validated)
- Sample: Serum or EDTA-plasma
- Assay time: appr. 3 h
- Required material not included: Adjustable micropipettes  
Microplate reader with measuring wavelength within 690-720 nm  
Clean containers for dilution of samples and chromogen

## Asparaginase Activity during Treatment with different Asparaginase Preparations using different Methods of Monitoring

ALL patient (♂/15 years) under Spectrila® treatment (5000 U/m<sup>2</sup>) during the first exposure (induction phase)




ALL patient (♂/2 years) under Asparaginase medac treatment (5000 U/m<sup>2</sup>) during the first exposure (induction phase)




# Development of Asparaginase Activity Levels

ALL patients treated with 1000 U/m<sup>2</sup> Oncaspar®<sup>6)</sup>




Plasma: asparaginase (ASP) activity levels\* (■) vs asparagine (ASN) levels\* (▲) during the first exposure (induction, phase IA) to Oncaspar® (given on days 12 and 27). \*mean ± SD



Plasma: asparaginase (ASP) activity levels\* (■) vs asparagine (ASN) levels\* (▲) during the second exposure (reinduction, protocols II or III) to PEG-ASP (considered to given conventionally given on day 1). \*mean ± SD.

- MAAT and AHA method deliver comparable results.
- Asparaginase levels  $\geq 100$  U/L correlate with complete asparagine depletion.
- MAAT can be used for various commercially available asparaginase preparations.<sup>3), 6), 8)</sup>

## Test Procedure



## Performance Characteristics

### Prevalence

102 sera obtained from blood donors were measured during the diagnostic evaluation. The asparaginase activity was generally < 30 U/L. In addition 52 sera obtained from children were investigated for asparaginase activity. The values of this population were generally < 30 U/L as well.

### Precision

**Intraassay Variation** determined in 4 reactive samples between 30-600 U/L, 22 fold measurement:

Sample	U/L	Mean OD	CV (%)
Calibrator	600.0	2.162	3.1
Control	96.4	0.530	1.2
Serum A	39.6	0.254	2.5
Serum B	310.3	1.380	1.5

**Interassay Variation** determined in 4 samples in 9 test runs, independent from each other:

Sample	Mean U/L	CV (%)
Control	91.5	4.1
Serum No. 1	38.2	4.5
Serum No. 2	420.3	8.0
Serum No. 3	125.5	6.0

## Recovery

A mean recovery of 93.6 % (SD = 9.9 %) was calculated adding 3 defined asparaginase activities (Asparaginase medac) each to 3 different sera.

## Dilution Linearity of reactive Sera

The dilution linearity was checked with highly reactive sera which were tested in 4-5 different dilutions (1:2 serial dilution) and showed a CV of 1-7 %.

	Dil. 1 U/L	Dil. 2 U/L	Dil. 3 U/L	Dil. 4 U/L	Dil. 5 U/L	Mean U/L	SD U/L	CV %
No. 1	1731	1789	1845	1897	-	1815	71	3.9
No. 2	440	446	465	476	-	457	17	3.6
No. 3	1569	1638	1665	1672	-	1636	47	2.9
No. 4	359	357	373	375	-	366	9	2.5
No. 5	373	380	407	404	-	391	17	4.4
No. 6	543	520	546	564	614	557	35	6.3
No. 7	940	945	972	979	-	959	19	2.0
No. 8	8876	8889	9248	8986	9125	9025	160	1.8
No. 9	3609	3775	4025	4273	-	3920	291	7.4
No. 10	1087	1068	1116	1153	-	1106	37	3.4

There was no hook effect observed for asparaginase activities up to 20000 U/L.

# MAAT

medac Asparaginase Activity Test

## Advantages of MAAT

MAAT is used as a clinical trial assay by several study groups in different trials worldwide (e.g. AIEOP-BFM ALL 2009; IntReALL 2010<sup>3</sup>); 58081 EORTC<sup>3, 8</sup>).

- Identification of patients with silent inactivation of asparaginase
- Suitable for measuring activity of various asparaginase preparations
- Optimizing of individual asparaginase dosing
- Drug monitoring at point of care
- Results available in 3 h
- High accuracy
- Easy-to-use: no washing or centrifugation steps

## References

- 1) Boos et al.: Monitoring of Asparaginase Activity and Asparagine Levels in Children on Different Asparaginase Preparations. *European Journal of Cancer*. 32A (9), 1544-1559 (1996).
- 2) Lanvers et al.: Analytical validation of a microplate reader-based method for the therapeutic drug monitoring of L-asparaginase in human serum. *Analytical Biochemistry*. 309, 117-126 (2002).
- 3) Leão and Costa: Therapeutic Drug Monitoring (TDM) of Asparaginase – the Porto Experience. Paper presented at The Children's Leukemia Group EORTC meeting; 2014 Oct 10-11; Porto, Portugal.
- 4) Mueller et al.: Pegylated Asparaginase (Oncaspar) in Children with ALL: Drug monitoring in reinduction according to the ALL/NHL-BFM protocols. *British Journal of Hematology*. 110, 379-384 (2000).
- 5) Riccardi et al.: L-Asparaginase pharmacokinetic and asparagine levels in cerebrospinal fluid of rhesus monkeys and humans. *Cancer Research*. 41, 4554-4558 (1981).
- 6) Rizzari et al.: A pharmacological study on pegylated Asparaginase used in front-line treatment of children with acute lymphoblastic leukemia. *Haematologica/The Hematology Journal*. 91(1), 24-31 (2006).
- 7) van der Sluis et al.: Consensus expert recommendations for identification and management of asparaginase hypersensitivity and silent inactivation. *Haematologica*. 101(3), 279-285 (2016).
- 8) Costa: Therapeutic Drug Monitoring (TDM) of Asparaginase – the Porto Experience. Poster session presented at the 10th Biennial Childhood Leukemia Symposium, 2016, Apr 25-26, Megaron, Athens, Greece.