


# EGFR Phospho Control Slides



### Intended Use

For In Vitro Diagnostic Use.

### Summary and Explanation

Epidermal Growth Factor Receptor is the receptor for epidermal growth factor. It is a member of the ErbB family receptors, a subfamily of four closely related receptor tyrosine kinases: EGFR (ErbB-1), HER-2 neu (ErbB-2), HER-3 (ErbB-3) and HER-4 (ErbB-4).

### Presentation

Five slides of EGFR Phospho positive tissues, each mounted on Hydrophilic Plus Slides, provided in a plastic mailer.

<i>Catalog No.</i>	<i>Quantity</i>
BSB-9167-CS	5 slides
BSB 6722	5 slides

**Storage** Store at 20-25°C

### Precautions

1. For professional users only. Results should be interpreted by a qualified medical professional.
2. Ensure proper handling procedures are used with this reagent.
3. Always wear personal protective equipment such as a laboratory coat, goggles, and gloves when handling reagents.
4. Dispose of unused solution with copious amounts of water.
5. Follow safety precautions of the heating device used for epitope retrieval (TintoRetriever Pressure Cooker or similar).
8. For additional safety information, refer to the Safety Data Sheet for this product.
9. For complete recommendations for handling biological specimens, please refer to the CDC document, "Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories" (see References in this document).

### Stability

**This product is stable up to the expiration date on the product label.**

Do not use the after expiration date listed on the package label.

### Slide Preparation

1. Deparaffinize, dehydrate and rehydrate tissues.
2. Subject tissues to enzymatic digestion with ImmunoDNA Digestor (BSB 0108 - BSB 0112) or similar product (Proteinase K) using the following protocol:
  - 2.1 Pre-warm the ImmunoDNA Digestor solution to room temperature.
  - 2.2 Incubate tissue sections with the ImmunoDNA Digestor solution for 10 minutes at room temperature or for 5 minutes at 37 °C.

### IHC Protocol

1. For manual staining, perform antibody incubation at ambient temperature. For automated staining methods, perform antibody incubation according to instrument manufacturer's instructions.
2. Wash slides with ImmunoDNA washer or DI water.
3. Continue IHC staining protocol. Wash slides between each step with ImmunoDNA washer solution.

### Abbreviated Immunohistochemical Protocol

Step	ImmunoDetector AP/HRP	PolyDetector AP/HRP	PolyDetector Plus HRP
Proteinase K	10 min	10 min	10 min
Peroxidase/AP Blocker	5 min.	5 min.	5 min
Primary Antibody	30-60 min.	30-60 min.	30-60 min.
1st Step Detection	10 min.	30-45 min.	15 min.
2nd Step Detection	10 min.	Not Applicable	15 min.
Substrate- Chromogen	5-10 min.	5-10 min.	5-10 min.
Counterstain / Coverslip	Varies	Varies	Varies

### Mounting Protocols

For detailed instructions using biodegradable permanent mounting media such as XyGreen PermaMunter (BSB 0169-0174) or organic solvent based resin such as PermaMunter (BSB 0094-0097), refer to PI0174 or PI0097.









### Product Limitations

Due to inherent variability present in immunohistochemical procedures (including fixation time of tissues, dilution factor of antibody, retrieval method utilized, and incubation time), optimal performance should be established through the use of positive and negative controls. Results should be interpreted by a qualified medical professional.

### References

1. Hunter T: Nature 1984, 311:414-416
2. Reynolds FH, Jr., et al.: Nature 1981, 292:259-262
3. Helin K, et al.: Oncogene 1991, 6:825-832
4. Christensen ME: Dan Med Bull 1998, 45:121-134
5. Kanematsu T, et al.: Oncol Res 2003, 13:289-298
6. Magkou C, et al.: Breast Cancer Res 2008, 10: R49
7. U.S. Department of Health and Human Services: Centers for Disease Control and Prevention. Guidelines for Safe <https://www.cdc.gov/mmwr/pdf/other/su6101.pdf>

### Symbol Key / Légende des symboles/Erläuterung der Symbole

	QAdvis EAR AB Ideon Science Park Scheelevägen 17 SE-223 70 Lund, Sweden	 Storage Temperature Limites de température Zulässiger Temperaturbereich	 Manufacturer Fabricant Hersteller	 Catalog Number Référence du catalogue Bestellnummer
	In Vitro Diagnostic Medical Device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	 Read Instructions for Use Lisez les instructions d'utilisation Gebrauchsanweisung beachten	 Expiration Date Utiliser jusque Verwendbar bis	 Lot Number Code du lot Chargenbezeichnung