

Rabbit anti-human AFP (Alpha-Fetoprotein) Monoclonal Antibody (Clone EP209)

REFERENCES AND PRESENTATIONS¹

- **ready-to-use (ml)**
MAD-000599QD-3
MAD-000599QD-7
MAD-000599QD-12
- **MD-Stainer presentations²**
MAD-000599QD-3/V
MAD-000599QD/V
- **concentrated**
MAD-000599Q - 1:50 recommended dilution

COMPOSITION

Anti-human AFP rabbit monoclonal antibody purified from serum and prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide

INTENDED USE : Immunohistochemistry (IHC) on paraffin embedded tissues. Not tested on frozen tissues or Western-Blotting

CLONE: EP209³

Ig ISOTYPE: Rabbit IgG

IMMUNOGEN: A synthetic peptide corresponding to residues of human AFP protein

SPECIES REACTIVITY: In vitro diagnostics in humans.

Not tested in other species

DESCRIPTION AND APPLICATIONS:

Alpha-fetoprotein (AFP) is the most abundant plasma protein found in the human fetus. It is thought to be the fetal form of serum albumin. AFP binds to copper, nickel, fatty acids and bilirubin and is found in monomeric, dimeric and trimeric forms.

Alpha-Fetoprotein (AFP) is synthesized by the cells of the embryonic yolk sac, fetal liver and fetal intestinal tract. AFP levels decrease soon after birth.

In abnormal tissues, expression of AFP has been demonstrated in hepatocellular carcinoma, hepatoid adenocarcinoma, germ cell tumors, particularly yolk sac tumor.

¹ These references are for presentation in vials of Low Density Polyethylene (LDPE) dropper. In case the products are used in automated stainings, a special reference is assigned as follows:

- / L: Cylindrical screw-cap vials (QD-3 / L, QD-7 / L, QD-12 / L).
- / N: Polygonal screw-cap vials (QD-3 / N, QD-7 / N, QD-12 / N).

For different presentations (references / volumes) please contact the supplier.

² For Technical specifications for MD-Stainer, please contact your distributor.

³ AFP clone EP209 is manufactured using Epitomics' RabMab® technology under U.S. Patent Nos. 5,675,063 and 7,402,409

The anti-AFP antibody may be useful for the identification of neoplastic liver diseases or yolk sac tumors.

IHC POSITIVE CONTROL: Foetal liver,

hepatocarcinoma or yolk sac tumour

VISUALIZATION: Cell cytoplasm

IHC RECOMMENDED PROCEDURE:

- 4µm thick section should be taken on charged slides; dry overnight at 60°C
- Deparaffinise, rehydrate and HIER (heat induced epitope retrieval) – boil tissue in the Pt Module using Master Diagnóstica EDTA buffer pH8⁴ for 20 min at 95°C. Upon completion rinse with 3-5 changes of distilled or deionised water followed by cooling at RT for 20 min
- Endogenous peroxidase block - Blocking for 10 minutes at room temperature using peroxidase solution (ref. MAD-021540Q-125)
- Primary antibody: incubate for 10 minutes [The antibody dilution (when concentrated) and protocol may vary depending on the specimen preparation and specific application. Optimal conditions should be determined by the individual laboratory]
- For detection use Master Polymer Plus Detection System (HRP) (DAB included; ref. MAD-000237QK)
- Counterstaining with haematoxylin and final mounting of the slide

STORAGE AND STABILITY:  up to 18 months;
 stored at 2-8°C. Do not freeze.

WARNINGS AND PRECAUTIONS:

1. Avoid contact of reagents with eyes and mucous membranes. If reagents come into contact with sensitive areas, wash with copious amounts of water.
2. This product is harmful if swallowed.
3. Consult local or state authorities with regard to recommended method of disposal.
4. Avoid microbial contamination of reagents.

SAFETY RECOMMENDATIONS

This product is intended for laboratory professional use only. The product is NOT intended to be used as a drug or for domestic purposes. The current version of the Safety Data Sheet for this product can be downloaded by searching the reference number at www.vitro.bio or can be requested at regulatory.md@vitro.bio.

⁴ Ref: MAD-004072R/D



BIBLIOGRAPHY

1. Nayak NC, Mital I. The dynamics of alpha-fetoprotein and albumin synthesis in human and rat liver during normal ontogeny. Am J Pathol 86: 359-74 (1977).
2. Smith CJP, Kelleher PC. Biochem Biophys Acta 650: 1-32. (1980).
3. Bosman FT, Giard RW, Nieuwenhuijen Kruseman AC, Knijnenburg G, Spaander PJ. Human chorionic gonadotrophin and alpha-fetoprotein in testicular germ cell tumours: a retrospective immunohistochemical study. Histopathology 4:673-84, (1980).
4. Jacobsen GK, Jacobsen M, Clausen PP. Distribution of tumor-associated antigens in the various histologic components of germ cell tumors of the testis. Am J Surg Pathol 5: 257-66, (1981).
5. Uchida T, Shimoda T, Miyata H, Shikata T, Iino S, Suzuki H, Oda T, Hirano K, Sugiura M. Immunoperoxidase study of alkaline phosphatase in testicular tumor. Cancer 48: 1455-62, (1981).
6. Richardson RL, Schoumacher RA, Fer MF, Hande KR, Forbes JT, Oldham RK, Greco FA. The unrecognized extragonadal germ cell cancer syndrome. Ann Int Med 94: 181-6, (1981).
7. Morinaga T, Sakai M, Wegmann TG, Tamaoki T. Primary structures of human alpha-fetoprotein and its mRNA. PNAS 80: 4604-8, (1983).
8. Imoto M, Nishimura D, Fukuda Y, Sugiyama K, Kumada T, Nakano S. Immunohistochemical detection of alpha-fetoprotein, carcinoembryonic antigen, and ferritin in formalin-paraffin sections from hepatocellular carcinoma. Am J Gastroenterol 80: 902-6, (1985).
9. Koelma IA, Nap M, Huitema S, Krom RA, Houthoff HJ. Hepatocellular carcinoma, adenoma, and focal nodular hyperplasia. Comparative histopathologic study with immunohistochemical parameters. Arch Pathol Lab Med 110: 1035-40, (1986).
10. Wick MR, Swanson PE, Manivel JC. Placental-like alkaline phosphatase reactivity in human tumors: an immunohistochemical study of 520 cases. Hum Pathol 18: 946-54, (1987).

