

**MASTER DIAGNÓSTICA, S.L.**

Avda. Conocimiento, 100

P.T. Ciencias de la Salud

18016 GRANADA

[www.masterdiagnostica.com](http://www.masterdiagnostica.com)

Fax: 958 27.14.34

Phone: 958 27.14.49

---

## MASTER GREASE for IHC

---

**APPLICATIONS:**

The MASTER GREASE marker prevents the reagents from spreading over the surface of the slide, creating a barrier that repels water. For this purpose, it is only necessary to make a circle around the tissue section. It is useful for immunohistochemical and immunofluorescence methods.

**DESCRIPTION:**

Hydrophobic linear marker with over 600 use applications.

**Reference: MAD-011520Q**

**RECOMMENDATIONS FOR USE:****1. FOR PARAFFIN EMBEDDED SECTIONS:**

- Deparaffinize and rehydrate the tissue sections.
- Wash the tissue section with PBS and drain the excess of liquid around the tissue.
- Make a circle on the slide around the tissue section with MASTER GREASE.
- Let the circle dry for 1 - 2 minutes at room temperature.
- Continue with the immunohistochemical staining process.

**2. FOR FROZEN SECTIONS:**

- Clean the slide.
- Cover the slide with an adhesive medium for frozen sections (ref. 011500Q).
- Cut and mount the frozen section on the slide.
- Wash the tissue section with PBS and drain the excess of liquid around the tissue.
- Make a circle on the slide around the tissue section with MASTER GREASE.
- Let the circle dry for 1 - 2 minutes at room temperature.
- Continue with the immunohistochemical staining process.

**NOTE:** The circle made with MASTER GREASE can be easily removed with xylene or a product for that purpose. It cannot be altered with water-soluble media and it does not interfere in microscopic observations of tissue sections.

**STORAGE:**

Store with the cap firmly closed and upside-down. Press the tip when the marker starts to get dry.

Keep refrigerated at 2°C - 8°C.

R: 60-36/37/38-48/20

S:53-26-45