














BioCLIA[®] Autoimmune Control Set

Jo-1

Jo-1 Assay Controls

Key to Symbols Used

	Catalog Number		Expiration Date
	For <i>In Vitro</i> Diagnostic Use		Lot Number
	Store between +2°C and +8°C		Consult Instruction for Use
	Manufacturer		Authorized Representative in European Union
	Control 1		Contains Sufficient for $< n >$ Tests
	Control 2		Chemical Risk Warning
			Biological Risk Warning

BioCLIA® Autoimmune Control Set,

Jo-1

Intended Use

The BioCLIA Autoimmune Control Set, Jo-1 is intended for the quality control purposes of the BioCLIA Jo-1 performed on the BioCLIA® 1200 and BioCLIA® 6500.

Catalog Numbers

MY00306 (2 X 1 mL)

MY00357 (4 X 1 mL)

Summary and Principles of the Procedure

Anti-nuclear antibody (ANA) is a class of auto-antibodies with different binding affinities specific to different nuclear antigens. Generally, ANAs include extractable nuclear antigen (ENA) antibodies and un-extractable nuclear antigen antibodies.¹ Determination of ANA has significant correlation to the diagnosis of Sharp syndrome (MCTD), systemic lupus erythematosus (SLE), sjogren's syndrome, progressive systemic sclerosis, polymyositis/dermatomyositis, overlap syndrome, and limited types of progressive systemic sclerosis (CREST syndrome).

Jo-1 is a 50 kDa cytoplasm phosphoprotein, the same component to histidyl-tRNA synthetase. It plays an important role in transporting free histidine to the corresponding tRNA and translating peptide from mRNA templates.²

Anti-Jo-1 antibodies are commonly seen in polymyositis (PM) patients' sera with the sensitivity about 40% ~ 50%, 25% in PM/DM (dermatomyositis) patients, less than 10% in DM patients, and not detected in other autoimmune disease patients. The sensitivity of anti-Jo-1 antibodies can be as high as 85% for the group who suffered PM and scleroderma, 25% for the ones with progressive systemic sclerosis (PSS)/PM and 60% for the PM ones with pulmonary fibrosis.^{3, 4}

Materials supplied

- **Jo-1 Control 1** A tube contains 1mL, ready to use reagent. Control contains human antibodies to Jo-1 in stabilizers and preservatives (low).

Preservatives: 0.0015% < Proclin 300 < 0.6%.

Jo-1	Control	L
------	---------	---

- **Jo-1 Control 2** A tube contains 1mL, ready to use reagent. Control contains human antibodies to Jo-1 in stabilizers and preservatives (high).

Preservatives: 0.0015% < Proclin 300 < 0.6%.

Jo-1	Control	H
------	---------	---

Target value and acceptable range for the controls are indicated on control information sheet in each kit.

Warnings and Precautions

The human derived material in this product was tested by FDA approved methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.⁵ Avoid contacting with skin and eyes. Do not empty into drains. Wear suitable protective clothing.

Precautions:



Human serum is added in the controls.



Proclin 300 is added in the controls at concentration between 0.0015% - 0.6%.

- The product is for *in vitro* diagnostic use only.
- Do not use any controls beyond their expiration dates. Do not mix controls from different lots unless specified.
- Instructions must be carefully followed for using and storing of controls. Any modification in procedure may interfere with the results. Controls and contaminated vials must be handled strictly following safety guidelines or rules of biological hazards to ensure the users' and environmental safety.
- Controls contain chemical and biological components. Avoid ingesting or splashing onto skin and mucous membrane. If direct contact with controls happens, rinse the contact surface with plenty of water immediately and see a doctor if necessary.

Storage Conditions

The kit is stable until the expiration date, if it is stored and handled as directed. Routine store the kit in refrigerator(2-8°C). Once a control tube is opened, it is good for a total of 15 times, no more than 2 hours per time when kept uncapped, onboard the instrument, after which the reagent must be discarded. Three freeze-thaw cycles before testing has no effect on the kit reagents.

Assay Procedure

Note that, for obtaining optimal performance, it is important to perform all routine maintenance procedures as defined in the BioCLIA® 1200 and BioCLIA® 6500 User Manual.

See the BioCLIA® 1200 and BioCLIA® 6500 User Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

The control procedure can be done before running the specimens each day. Users also can adjust the control procedure period according to their own lab frequency.

Limitations

This product is designed as controls for monitoring

the performance of the BioCLIA Jo-1. These controls are subjected to the limitations of the assay system. Deviations may indicate possible problems with one or more components in the test system.

References

1. Tan EM. Autoantibodies to nuclear antigens (ANA): their immunobiology and medicine. *Advances in Immunology* 1982;33:167-240.
2. Damoiseaux J, Tervaert JC. From ANA to ENA: how to proceed? *Autoimmunity reviews* 2006;5:10-17.
3. Autoantibodies in the diagnosis of systemic rheumatic diseases. *Seminars in arthritis and rheumatism*; 1995. Elsevier.
4. Lyons R, Narain S, Nichols C, Satoh M, Reeves WH. Effective use of autoantibody tests in the diagnosis of systemic autoimmune disease. *Annals of the New York Academy of Sciences* 2005;1050:217-28.
5. Richmond JY, McKinney RW. Biosafety in microbiological and biomedical laboratories: U.S.GPO. 1999.



HOB Biotech Group Co., Ltd

C6 Building, No. 218 Xinghu Road, Suzhou Industrial Park,

Suzhou, Jiangsu, 215123, China

REGISTRANT/MANUFACTURE: HOB Biotech Group Co., Ltd

ADDRESS/LOCATION:

C6 Building, No. 218 Xinghu Road, Suzhou Industrial Park, Suzhou, Jiangsu, 215123 China

CONTACT INFORMATION: TEL (+86)512-69561996
Fax (+86)512-62956652

WEBSITE: www.hob-biotech.com

CUSTOMER SERVICE: HOB Biotech Group Co., Ltd

CUSTOMER SERVICE CONTACT: TEL (+86)4008601202



EUROPE REPRESENTATIVE: Emergo Europe

ADDRESS/LOCATION:

Prinsessegracht 20, 2514 AP The Hague, The Netherlands

Technical Assistance

For technical assistance, contact your National Distributor.

17th April 2019

Revision 9