














BioCLIA[®] Autoimmune Control Set

PCNA

PCNA Assay Controls

Key to Symbols Used			
	Catalog Number		Expiration Date
	For <i>In Vitro</i> Diagnostic Use		Lot Number
	Store between +2°C and +8°C		Consult Instruction for Use
	Manufacturer		Authorized Representative in European Union
	Control 1		Contains Sufficient for n Tests
	Control 2		Chemical Risk Warning
			Biological Risk Warning

BioCLIA® Autoimmune Control Set, PCNA

Intended Use

The BioCLIA Autoimmune Control Set, PCNA is intended for the quality control purposes of the BioCLIA PCNA performed on the BioCLIA®1200 and BioCLIA® 6500.

Catalog Numbers

MY00310 (2 X 1 mL)

My00361 (4 X 1 mL)

Summary and Explanation

Anti-nuclear antibodies (ANAs) are a class of auto-antibodies with different binding affinities specific to different nuclear antigens. Generally, ANAs include extractable nuclear antigen (ENA) antibodies and un-extractable nuclear antigen antibodies.¹ Determination of ANAs has significant correlation to the diagnosis of Sharp Syndrome, Systemic Lupus Erythematosus (SLE), Sjogren's Syndrome, progressive systemic sclerosis, polymyositis/dermatomyositis, overlap syndrome, and limited types of progressive systemic sclerosis (CREST syndrome).

Proliferating cell nuclear antigen (PCNA) is a DNA clamp that acts as a processivity factor for DNA polymerase δ in eukaryotic cells and is essential for replication. Antibodies against proliferating cell nuclear antigen (PCNA, a kind of ANA marker) or monoclonal antibody termed Ki-67 can be used for grading of different neoplasms, e.g. astrocytoma. They can be of diagnostic and prognostic value. Imaging of the nuclear distribution of PCNA (via antibody labeling) can be used to distinguish between early, mid and late S phase of the cell cycle.²⁻⁴

Materials supplied

- **PCNA Control 1** A tube contains 1mL, ready to use reagent. Control contains human antibodies to PCNA in stabilizers and preservatives (Low).

PCNA | Control | L

Preservatives: 0.0015% < Proclin 300 < 0.6%.

- **PCNA Control 2** A tube contains 1mL, ready to use reagent. Control contains human antibodies to PCNA in stabilizers and preservatives (High).

PCNA | Control | H

Preservatives: 0.0015% < Proclin 300 < 0.6%.

Target value and acceptable range for the controls are indicated on control information sheet in each kit.

Warnings and Precautions

The human derived material in this product was tested by FDA approved methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.⁵ Avoid contacting with skin and

eyes. Do not empty into drains. Wear suitable protective clothing.

Precautions:



Human serum is added in the controls.



Proclin 300 is added in the controls at concentration between 0.0015% - 0.6%.

- The product is for *in vitro* diagnostic use only.
- Do not use any controls beyond their expiration dates. Do not mix controls from different lots unless specified.
- Instructions must be carefully followed for using and storing of controls. Any modification in procedure may interfere with the results. Controls and contaminated vials must be handled strictly following safety guidelines or rules of biological hazards to ensure the users' and environmental safety.
- Controls contain chemical and biological components. Avoid ingesting or splashing onto skin and mucous membrane. If direct contact with controls happens, rinse the contact surface with plenty of water immediately and see a doctor if necessary.

Storage Conditions

The kit is stable until the expiration date, if it is stored and handled as directed. Routine store the kit in refrigerator(2-8°C). Once a control tube is opened, it is good for a total of 15 times, no more than 2 hours per time when kept uncapped, onboard the instrument, after which the reagent must be discarded. Three freeze-thaw cycles before testing has no effect on the kit reagents.

Assay Procedure

Note that, for obtaining optimal performance, it is important to perform all routine maintenance procedures as defined in the BioCLIA®1200 and BioCLIA® 6500 User Manual.

See the BioCLIA®1200 and BioCLIA® 6500 User Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

The control procedure can be done before running the specimens each day. Users also can adjust the control procedure period according to their own lab frequency.

Limitations

This product is designed as controls for monitoring the performance of the BioCLIA PCNA. These controls are subjected to the limitations of the assay system. Deviations may indicate possible problems with one or more components in the test system.

References

1. Tan EM. Autoantibodies to nuclear antigens (ANA): their immunobiology and medicine. *Advances in Immunology* 1982;33:167-240.
2. Schönenberger F, Deutzmann A, Ferrandomay E, Merhof D. Discrimination of cell cycle phases in PCNA-immunolabeled cells. *BMC Bioinformatics* 2015;16:1-10.
3. Herce HD, Rajan M, Lättig-Tünnemann G, Fillies M, Cardoso MC. A novel cell permeable DNA replication and repair marker. *Nucleus* 2014;5:590-600.
4. Wang S-C. PCNA: a silent housekeeper or a potential therapeutic target? *Trends in pharmacological sciences* 2014;35:178-86.
5. Richmond JY, Mckinney RW. Biosafety in microbiological and biomedical laboratories: U.S.GPO. 1999.



HOB Biotech Group Co., Ltd
C6 Building, No. 218 Xinghu Road, Suzhou Industrial Park,
Suzhou, Jiangsu, 215123, China

REGISTRANT/MANUFACTURE: HOB Biotech Group Co., Ltd

ADDRESS/LOCATION:

C6 Building, No. 218 Xinghu Road, Suzhou Industrial Park, Suzhou, Jiangsu, 215123 China

CONTACT INFORMATION: TEL (+86)512-69561996
Fax (+86)512-62956652

WEBSITE: www.hob-biotech.com

CUSTOMER SERVICE: HOB Biotech Group Co., Ltd

CUSTOMER SERVICE CONTACT: TEL (+86)4008601202



EUROPE REPRESENTATIVE: Emergo Europe

ADDRESS/LOCATION:

Prinsessegracht 20, 2514 AP The Hague, The Netherlands

Technical Assistance

For technical assistance, contact your National Distributor.

17th April 2019

Revision 9