


## CD2

### Control Slides



#### Intended Use

For In Vitro Diagnostic Use.

#### Summary and Explanation

CD2 is a cell-adhesion molecule found on the surface of T-cells and natural killer cells. It has also been called T-cell surface antigen T11/Leu-5, LFA-2, LFA-3 receptor, erythrocyte receptor and rosette receptor. Due to its structural characteristics, CD2 is a member of the immunoglobulin superfamily; it possesses two immunoglobulin-like domains in its extracellular portion. It interacts with other adhesion molecules, such as lymphocyte function-associated antigen-3 (LFA-3/CD58) in humans, or CD48 in rodents, which are expressed on the surfaces of other cells. In addition to its adhesive properties, CD2 also acts as a co-stimulatory molecule on T and NK cells.

CD2 is a surface antigen of the human T-lymphocyte lineage that is expressed on all peripheral blood T-cells. It is one of the earliest T-cell markers, being present on more than 95% of thymocytes; it is also found on some natural killer cells but not on B-lymphocytes. CD2 is implicated in the triggering of T-cells; the cytoplasmic domain is implicated in the signaling function. It is useful for the identification of Lymphomas and Leukemias of T-cell origin. As with other pan-T cell antigens, CD2 may be aberrantly deleted in some neoplastic T-cell populations, especially Peripheral T-cell Lymphomas.

#### Presentation

Five slides of CD2 positive tissues, each mounted on Hydrophilic Plus Slides, provided in a plastic mailer.

| <i>Catalog No.</i> | <i>Quantity</i> |
|--------------------|-----------------|
| BSB-9077-CS        | 5 slides        |
| BSB 6211           | 5 slides        |

**Storage** Store at 20-25°C

#### Precautions

1. For professional users only. Results should be interpreted by a qualified medical professional.
2. Ensure proper handling procedures are used with this reagent.
3. Always wear personal protective equipment such as a laboratory coat, goggles, and gloves when handling reagents.
4. Dispose of unused solution with copious amounts of water.
5. Follow safety precautions of the heating device used for epitope retrieval (TintoRetriever Pressure Cooker or similar).
8. For additional safety information, refer to Safety Data Sheet for this product.
9. For complete recommendations for handling biological specimens, please refer to the CDC document, "Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories" (see References in this document).

#### Stability

**This product is stable up to the expiration date on the product label.**

Do not use after expiration date listed on package label.

#### IHC Protocol

1. Subject tissues to heat induced epitope retrieval (HIER) using a suitable retrieval solution such as ImmunoDNA Retriever with Citrate (BSB 0020-BSB 0023) or EDTA (BSB 0030-BSB 0033).

2. Any of three heating methods may be used:

##### **a. TintoRetriever Pressure Cooker or Equivalent**

Place tissues/slides in a staining dish or coplin jar containing the ImmunoDNA Retriever with Citrate or EDTA and place on trivet in the pressure cooker. Add 1-2 inches of distilled water to the pressure cooker and turn heat to high. Incubate for 15 minutes. Open and immediately transfer slides to room temperature.

##### **b. TintoRetriever PT Module or Water Bath Method**

Place tissues/slides in a pre-warmed staining dish or coplin jar containing the ImmunoDNA Retriever with Citrate or EDTA at 95°-99° C. Incubate for 30-60 minutes.

##### **c. Conventional Steamer Method**

Place tissues/slides in a pre-warmed staining dish or coplin jar containing the ImmunoDNA Retriever with Citrate or EDTA in a steamer, cover and steam for 30-60 minutes.

3. After heat treatment, transfer slides in ImmunoDNA Retriever with Citrate or EDTA to room temperature and let stand for 15-20 minutes.
4. For manual staining, perform antibody incubation at ambient temperature. For automated staining methods, perform antibody incubation according to instrument manufacturer's instructions.
5. Wash slides with ImmunoDNA washer or DI water.
6. Continue IHC staining protocol. Wash slides between each step with ImmunoDNA washer solution.

### Abbreviated Immunohistochemical Protocol

| Step                     | ImmunoDetector<br>AP/HRP | PolyDetector<br>AP/HRP | PolyDetector<br>Plus HRP |
|--------------------------|--------------------------|------------------------|--------------------------|
| Peroxidase/AP Blocker    | 5 min.                   | 5 min.                 | 5 min.                   |
| Primary Antibody         | 30-60 min.               | 30-60 min.             | 30-60 min.               |
| 1st Step Detection       | 10 min.                  | 30-45 min.             | 15 min.                  |
| 2nd Step Detection       | 10 min.                  | Not Applicable         | 15 min.                  |
| Substrate- Chromogen     | 5-10 min.                | 5-10 min.              | 5-10 min.                |
| Counterstain / Coverslip | Varies                   | Varies                 | Varies                   |

### Abbreviated IF Protocol

| Step   | Incubation Time   |
|--|-------------------|
| Rinse slides in IF wash buffer                 | 5 minutes         |
| Drain and wipe excess IF wash buffer off slide |                   |
| Conduct remaining steps in the dark            |                   |
| Apply Antibody                                 | 30-60 minutes     |
| Rinse with 3 changes of IF wash buffer         | 3x15 minutes each |
| Coverslip with IF mounting medium              |                   |

### Mounting Protocols

For detailed instructions using biodegradable permanent mounting media such as XyGreen PermaMunter (BSB 0169-0174) or organic solvent based resin such as PermaMunter (BSB 0094-0097), refer to PI0174 or PI0097.

### Product Limitations

Due to inherent variability present in immunohistochemical procedures (including fixation time of tissues, dilution factor of antibody, retrieval method utilized, and incubation time), optimal performance should be established through the use of positive and negative controls. Results should be interpreted by a qualified medical professional.

### References

- Hyjek E, Chadburn A, Liu YF, Cesarman E, Knowles DM. BCL-6 protein is expressed in precursor T-cell lymphoblastic lymphoma and in prenatal and postnatal thymus. Blood 2001;97:270-76.
- Went P, Agostinelli C, Gallamini A, Piccaluga PP, Ascani S, Sabattini E, et al. Marker expression in peripheral T-cell lymphoma: A proposed clinical-pathological prognostic score. J Clin Oncol 2006;24:2472-79.
- D'Amore SEG, Menin A, Bonoldi E, Bevilacqua P, Cazzavilan S, Donofrio V, et al. Anaplastic large cell lymphomas: A study of 75 pediatric patients. Pediatr Dev Pathol 2007;10:181-91.
- Leong AS-Y, Cooper K and Leong FJW-M. CD2. Manual of diagnostic antibodies for immunohistology. London: Greenwich Medical Media; 2003. p. 61-62.
- Moingeon P, Chang HC, Sayre PH, Clayton LK, Alcover A, Gardner P, et al. The structural biology of CD2. Immunol Rev 1989;111:111-44.
- U.S. Department of Health and Human Services: Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. Supplement / Vol. 61, January 6, 2012. <https://www.cdc.gov/mmwr/pdf/other/su6101.pdf>

### Symbol Key / Légende des symboles/Erläuterung der Symbole

|                      |  |   |   |            |   |
|----------------------|--|---|---|------------|---|
| <b>EC</b> <b>REP</b> | QAdvis EAR AB<br>Ideon Science Park<br>Scheelevägen 17<br>SE-223 70 Lund, Sweden                         |  Storage Temperature<br>Limites de température<br>Zulässiger Temperaturbereich                           |  Manufacturer<br>Fabricant<br>Hersteller               | <b>REF</b> | Catalog Number<br>Référence du catalogue<br>Bestellnummer |
| <b>IVD</b>           | In Vitro Diagnostic Medical Device<br>Dispositif médical de diagnostic in vitro<br>In-Vitro-Diagnostikum |  Read Instructions for Use<br>Consulter les instructions<br>d'utilisation<br>Gebrauchsanweisung beachten |  Expiration Date<br>Utiliser jusqu'à<br>Verwendbar bis | <b>LOT</b> | Lot Number<br>Code du lot<br>Chargenbezeichnung           |