


p14 ARF Control Slides



Intended Use

For In Vitro Diagnostic Use.

Summary and Explanation

p14 ARF (also called ARF tumor suppressor, ARF, p14ARF) encoded by the p16 tumor suppressor gene is an alternate reading frame protein product of the CDKN2A locus (i.e. INK4a/ARF locus). p14 ARF accumulates mainly in the nucleolus where it forms stable complexes with NPM or MDM2. These interactions allow p14 ARF to act as a tumor suppressor by inhibiting ribosome biogenesis or initiating p53-dependent cell cycle arrest and apoptosis, respectively. Both p16INK4a and p14 ARF are involved in cell cycle regulation. p14 ARF inhibits MDM2, thus promoting p53, which promotes p21 activation, which then binds and inactivates certain cyclin-CDK complexes, which would otherwise promote transcription of genes that would carry the cell through the G1/S checkpoint of the cell cycle. Loss of p14 ARF by a homozygous mutation in the CDKN2A (INK4A) gene will lead to elevated levels in MDM2 and, therefore, loss of p53 function and cell cycle control. p14 ARF, has been reported to be associated with the clinicopathological features of different cancers. Very commonly, cancer is associated with a loss of function of INK4a, ARF, Rb, or p53. Without ARF, MDM2 can inappropriately inhibit p53, leading to increased cell survival. The INK4a/ARF locus is found to be deleted or silenced in many kinds of tumors. It has been found that 41% breast carcinomas have p14 ARF defects and in a separate study, 32% of colorectal adenomas were found to have p14 ARF inactivation due to hypermethylation of the promoter. Homozygous deletions and other mutations of CDK2NA (ARF) have been found to be associated with Glioblastoma. p14 ARF expression has been found to be significantly associated with the risk of lung cancer.

Presentation

Five slides of p14 ARF positive tissues, each mounted on Hydrophilic Plus Slides, provided in a plastic mailer.

<i>Catalog No.</i>	<i>Quantity</i>
BSB-9320-CS	5 slides
BSB 3601	5 slides

Storage Store at 20-25°C

Precautions

1. For professional users only. Results should be interpreted by a qualified medical professional.
2. Ensure proper handling procedures are used with this reagent.
3. Always wear personal protective equipment such as a laboratory coat, goggles, and gloves when handling reagents.
4. Dispose of unused solution with copious amounts of water.
5. Follow safety precautions of the heating device used for epitope retrieval (TintoRetriever Pressure Cooker or similar).
8. For additional safety information, refer to Safety Data Sheet for this product.
9. For complete recommendations for handling biological specimens, please refer to the CDC document, "Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories" (see References in this document).

Stability

This product is stable up to the expiration date on the product label.

Do not use after expiration date listed on package label.

IHC Protocol

1. Subject tissues to heat induced epitope retrieval (HIER) using a suitable retrieval solution such as ImmunoDNA Retriever with Citrate (BSB 0020-BSB 0023) or EDTA (BSB 0030-BSB 0033).

2. Any of three heating methods may be used:

a. TintoRetriever Pressure Cooker or Equivalent

Place tissues/slides in a staining dish or coplin jar containing the ImmunoDNA Retriever with Citrate or EDTA and place on trivet in the pressure cooker. Add 1-2 inches of distilled water to the pressure cooker and turn heat to high. Incubate for 15 minutes. Open and immediately transfer slides to room temperature.

b. TintoRetriever PT Module or Water Bath Method

Place tissues/slides in a pre-warmed staining dish or coplin jar containing the ImmunoDNA Retriever with Citrate or EDTA at 95°-99° C. Incubate for 30-60 minutes.

c. Conventional Steamer Method

Place tissues/slides in a pre-warmed staining dish or coplin jar containing the ImmunoDNA Retriever with Citrate or EDTA in a steamer, cover and steam for 30-60 minutes.

3. After heat treatment, transfer slides in ImmunoDNA Retriever with Citrate or EDTA to room temperature and let stand for 15-20 minutes.
4. For manual staining, perform antibody incubation at ambient temperature. For automated staining methods, perform antibody incubation according to instrument manufacturer's instructions.
5. Wash slides with ImmunoDNA washer or DI water.
6. Continue IHC staining protocol. Wash slides between each step with ImmunoDNA washer solution.

Abbreviated Immunohistochemical Protocol

Step	ImmunoDetector AP/HRP	PolyDetector AP/HRP	PolyDetector Plus HRP
Peroxidase/AP Blocker	5 min.	5 min.	5 min
Primary Antibody	30-60 min.	30-60 min.	30-60 min.
1st Step Detection	10 min.	30-45 min.	15 min.
2nd Step Detection	10 min.	Not Applicable	15 min.
Substrate- Chromogen	5-10 min.	5-10 min.	5-10 min.
Counterstain / Coverslip	Varies	Varies	Varies

Abbreviated IF Protocol

Step	Incubation Time
Rinse slides in IF wash buffer	5 minutes
Drain and wipe excess IF wash buffer off slide	
Conduct remaining steps in the dark	
Apply Antibody	30-60 minutes
Rinse with 3 changes of IF wash buffer	3x15 minutes each
Coverslip with IF mounting medium	

Mounting Protocols

For detailed instructions using biodegradable permanent mounting media such as XyGreen PermaMunter (BSB 0169-0174) or organic solvent based resin such as PermaMunter (BSB 0094-0097), refer to PI0174 or PI0097.

Product Limitations

Due to inherent variability present in immunohistochemical procedures (including fixation time of tissues, dilution factor of antibody, retrieval method utilized, and incubation time), optimal performance should be established through the use of positive and negative controls. Results should be interpreted by a qualified medical professional.

References

- Sherr CJ. "Divorcing ARF and p53: an unsettled case". Nat. Rev. Cancer. 2006; 6 (9): 663-73.
- Abida WM, Gu W (January 2008). "p53-Dependent and p53-independent activation of autophagy by ARF". Cancer Res. 2008; 68 (2): 352-7
- Sherr CJ (May 2006). "Autophagy by ARF: a short story". Mol. Cell. 2006; 22 (4): 436-7.
- Lowe SW, Sherr CJ. "Tumor suppression by Ink4a-Arf: progress and puzzles". Curr. Opin. Genet. Dev. 2003; 13 (1): 77-83.
- Yi Y, Shepard A, Kittrell F, Mulac-Jericevic B, Medina D, Said TK (May 2004). "p19ARF Determines the Balance between Normal Cell Proliferation Rate and Apoptosis during Mammary Gland Development". Mol. Biol. Cell. 2004; 15 (5): 2302-11.
- Cancer Genome Atlas Research, Network. "Comprehensive genomic characterization defines human glioblastoma genes and core pathways". Nature. 2008; 455(7216): 1061-8.
- Fang Wang, et al. Clinicopathological significance of p14ARF expression in lung cancer: a meta-analysis. OncoTargets and Therapy 2017;10 2491-2499
- U.S. Department of Health and Human Services: Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. Supplement / Vol. 61, January 6, 2012
<https://www.cdc.gov/mmwr/pdf/other/su6101.pdf>

Symbol Key / Légende des symboles/Erläuterung der Symbole

EC REP	QAdvis EAR AB Ideon Science Park Scheelevägen 17 SE-223 70 Lund, Sweden	 Storage Temperature Limites de température Zulässiger Temperaturbereich	 Manufacturer Fabricant Hersteller	REF	Catalog Number Référence du catalogue Bestellnummer
IVD	In Vitro Diagnostic Medical Device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	 Read Instructions for Use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Expiration Date Utiliser jusqu'à Verwendbar bis	LOT	Lot Number Code du lot Chargenbezeichnung