

# Peroxidase-Labeled Protein A from *S. aureus*



**Catalog No.** 14-50-00  
**Size** 0.5 mg

## Description

Chromatographically purified Protein A obtained from *Staphylococcus aureus* and labeled with horseradish peroxidase.

## FORM/STORAGE

Lyophilized. Store at 2 – 8°C until rehydrated. Stable for a minimum of one year at 2 – 8°C.

## STABILIZER AND PRESERVATIVE

Bovine serum albumin (BSA) added as protein stabilizer. No preservatives added. Additional biological protection may be provided with 0.01% thimerosal. DO NOT USE SODIUM AZIDE. Non-sterile.

## E/P RATIO

Molar enzyme/antibody protein ratio = 2:1

## SPECIFICITY/CROSS-REACTIVITY

Peroxidase-labeled Protein A adheres to the IgG of most animals<sup>1</sup> (Table 1). The affinity of Protein A to immunoglobulins increases with increasing pH to pH 8.5 and with increasing ionic strength to 0.3M. Low concentrations (0.05%) of Tween 20 markedly reduce the affinity of Protein A for IgG. To obtain optimal performance of this reagent, the user must adjust the pH, conductivity and detergent level to the needs of his experiment. KPL tests each lot of Protein A for reactivity to human, rabbit, mouse and swine IgG only.

TABLE 1

| SPECIES    | CLASSES/SUBCLASSES THAT PROTEIN A BINDS  |
|------------|--|
| Man        | IgG <sub>1</sub> , IgG <sub>2</sub> , IgG <sub>4</sub> and IgA <sub>2</sub>                  |
| Rabbit     | IgG and some IgM   |
| Mouse      | IgG <sub>2a</sub> , IgG <sub>2b</sub> , IgG <sub>3</sub> and some IgG <sub>1</sub>           |
| Sheep      | IgG <sub>2</sub>   |
| Goat       | IgG <sub>2</sub>   |
| Rat        | IgG <sub>1</sub> and IgG <sub>2c</sub>   |
| Dog        | IgG <sub>a</sub> , IgG <sub>b</sub> , IgG <sub>c</sub> , IgG <sub>d</sub> , IgA and some IgM |
| Guinea Pig | IgG <sub>1</sub> and IgG <sub>2</sub>  |
| Pig        | IgG  |
| Cow        | IgG <sub>2</sub>   |

## REHYDRATION AND STORAGE

### Rehydration:

Rehydrate with 1 mL of reagent quality water. Dilute to desired concentration with PBS or other buffer.

## Storage:

This product may be stored for up to 1 week refrigerated; thereafter, it should be stored frozen. When frozen, product is stable for a minimum of 1 year.

NOTE: Variations in temperature due to freeze cycles may cause loss of activity when rehydrated product is stored frozen in aliquots less than 50 µL.

## SUGGESTED WORKING CONCENTRATION

Different assay conditions require that serial dilutions of all reagents be performed to determine optimal working concentrations. Prepare working dilution immediately before use. Storage at working dilution may result in enzyme inactivation and performance loss. Suggested starting concentrations are as follows:

Microwell Plate Immunoassays: 1/500 (1.0 µg/mL)  
Membrane Immunoassays: 1/2500 (0.2 µg/mL)  
Histo/Cytochemical Procedures: 1/100 (5.0 µg/mL)

## CAUTION

Horseradish peroxidase is inactivated in the presence of hydrogen peroxide by reacting irreversibly with certain pollutants common in laboratory water supplies. If this product fails to perform as expected, check water supply for enzyme inactivation.

## PRODUCT SAFETY AND HANDLING

This product is considered non-hazardous as defined by The Hazard Communications Standard (29 CFR 1910.1200). Avoid contact with skin and eyes. In case of contact or spillage, clean with copious amounts of water. Disposal via sanitary sewer.

## RELATED PRODUCTS

10% BSA Diluent/Blocking Solution Catalog No. 50-61-00  
Coating Solution Concentrate Catalog No. 50-84-00  
Milk Diluent/Blocking Solution Catalog No. 50-82-01  
Wash Solution Concentrate Catalog No. 50-63-00

## Suggested KPL Peroxidase Substrates

### ELISA:

ABTS Microwell Substrate Catalog No. 50-66-00  
SureBlue Microwell Substrate Catalog No. 52-00-01

### Immunoblotting:

4CN Substrate System Catalog No. 50-73-00  
TMB Membrane Substrate Catalog No. 50-77-03

**Immunohistochemistry:**

|                     |                      |
|---------------------|----------------------|
| HistoMark® TrueBlue | Catalog No. 54-78-00 |
| HistoMark® ORANGE   | Catalog No. 54-74-00 |
| HistoMark® BLACK    | Catalog No. 54-75-00 |

See the KPL catalog for a wide selection of antibodies, substrate, protein and nucleic acid detection kits and immunohistochemistry reagents.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However, said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product. No claims beyond replacement of unacceptable material or refund of purchase price shall be allowed.