

PhosphaGLO Reserve™ AP Substrate



<u>Catalog No.</u>	<u>Size</u>
55-60-02	100 mL
55-60-01	30 mL

DESCRIPTION

PhosphaGLO Reserve AP Substrate contains a dioxetane-based chemiluminescent substrate designed for use with phosphatase-labeled (AP) reporter molecules. PhosphaGLO Reserve offers great improvements in the way of signal intensity. These products are specifically designed for the detection of proteins that are (1) difficult to detect because they are in such low quantities or (2) are from samples that are precious.

PhosphaGLO Reserve AP Substrate is provided as a stable one-component solution. The product provides rapid and accurate identification of proteins that are of low abundance and potentially limited availability. Given the increased sensitivity, less target may be required on Western blots and ELISAs.

Two sizes are available. Results can be obtained on X-ray film or a chemiluminescent imager to provide a permanent record. A luminometer should be used when performing ELISA.

NOTE: This product may arrive frozen.

CONTENTS

Catalog No. 55-60-01 contains:

1 x 30 mL PhosphaGLO Reserve AP Substrate

Catalog No. 55-60-02 contains:

1 x 100 mL PhosphaGLO Reserve AP Substrate

STORAGE/STABILITY

Store this at 2 – 8°C. If the product arrives frozen, allow it to thaw prior to use. PhosphaGLO Reserve Solution should remain stored in its original container and protected from light.

The product is stable for a minimum of two years from date of receipt when stored under proper conditions.

PRODUCT SAFETY AND HANDLING

See MSDS (Material Safety Data Sheet) for this product.

RELATED PRODUCTS

<u>Product</u>	<u>Catalog No.</u>
PhosphaGLO AP Substrate, 30 mL	55-60-04
PhosphaGLO AP Substrate, 100 mL	55-60-03
Phosphatase Assay Buffer (10X), 200 mL	50-63-14
Detector™ Block, 240 mL	71-83-00
Milk Diluent/Blocking Solution Concentrate Kit, 200 mL	50-82-01
10% BSA Diluent/Blocking Solution Kit, 200 mL	50-61-00
10% BSA Diluent/Blocking Solution, 1 L	50-61-10
Wash Solution Concentrate (20X), 800 mL	50-63-00
Coating Solution Concentrate, 50 mL	50-84-00

PHOSPHAGLO RESERVE AP SUBSTRATE USER'S GUIDE

- ⇒ PhosphaGLO Reserve can be used with nitrocellulose and PVDF membranes. For best results, nitrocellulose is recommended.
- ⇒ PhosphaGLO Reserve should be protected from light and warmed to room temperature prior to use.
- ⇒ PhosphaGLO Reserve is an extremely sensitive substrate. Insufficient washing of membranes or contamination of substrate with AP will result in non-specific background.
- ⇒ Because of PhosphaGLO Reserve's super sensitivity, the AP conjugate should be titrated to give the best results.
- ⇒ Do not allow PhosphaGLO Reserve to contact the film. If this occurs, PhosphaGLO Reserve solution will cause dark spots to appear on the film.
- ⇒ PhosphaGLO Reserve emits light over the course of 5-7 days. Because of its high light intensity, images may be captured over the course of many different intervals. Initially, perform a 1 minute and a 10 minute exposure.

APPLICATIONS

PhosphaGLO Reserve AP Substrate has been optimized for Western blotting and dot blotting applications. It is also suitable for use in microwell applications such as ELISA. The following is a recommended procedure for Western blot detection.

WESTERN BLOT DETECTION

Suggested Reagents/Equipment Not Included

1. Primary antibody
2. AP-labeled secondary antibody
3. Nitrocellulose or PVDF membrane
4. Blocking Solution (See RELATED PRODUCTS)
5. X-ray film (double emulsion) or CCD Imager
6. Platform shaker or rocker
7. Developing chemicals/equipment
8. Incubation trays or tubes
9. Wash solution (see RELATED PRODUCTS)
10. Assay buffer (see RELATED PRODUCTS)

CONJUGATE OPTIMIZATION PRIOR TO DETECTION

- ⇒ Before beginning the assay, the optimal conjugate dilution should be determined for the assay. The use of highly sensitive chemiluminescent substrates on Western blots can cause high background if the conjugate concentration is not optimized. Each lot of conjugate will need optimization as slight differences in activity can result in major differences in background.
- ⇒ Recommended conjugate dilutions should be tested in a range from 1/10,000 to 1/100,000 of a 0.1 mg/mL stock.

WESTERN BLOT DETECTION AT A GLANCE

Total time: 4 hours

Polyacrylamide Gel Electrophoresis



Immobilize Protein on Membrane



Block Membrane

1 hour or overnight



Incubate with Primary Antibody

1 hour



Wash Membrane

3 x 5 minutes

1 x 10 minutes



Incubate with Conjugate

1 hour



Wash Membrane

3 x 5 minutes

1 x 10 minutes



Rinse with Assay buffer

2x2 minutes



Incubate with PhosphaGLO Reserve Substrate

1 minute



Expose to Film

10 seconds - 10 minutes

STEPS	CRITICAL POINTS	STEPS	CRITICAL POINTS
<p>1. Block the membrane by immersing in block solution (1X Detector Block is recommended) using a minimum of 0.2 mL/cm² of membrane. Block at room temperature for 1 hour with gentle rocking, or stationary at 2 - 8°C overnight.</p>	<p><i>Example: for a 10 x 10 cm blot, use 20 mL of block solution. Make sure to use a container of proper size that allows the block solution to float freely over the membrane.</i></p>	<p>7. After the conjugate incubation, wash as described in step 3.</p>	
<p>2. Incubate membrane with primary antibody or serum sample for at least 1 hour. This antibody should be added directly to the Block Solution that was used for blocking (Step 1).</p>	<p><i>Test serial dilutions through a dot blot to determine the optimal working dilution or use the recommended concentration determined by the primary antibody supplier.</i></p>	<p>8. Pour off remaining wash buffer from the blot and place the membrane on a sheet protector or a dry tray.</p>	
<p>3. Wash the membrane in a generous amount of 1X Wash Solution (at least 25 mL for a 100 cm² membrane). Wash membrane 3 times for 5 minutes each, followed by one 10-minute wash.</p>	<p><i>1X TBS-TWEEN™ may also be used.</i></p>	<p>9. Rinse membrane 2 X 2 minutes with 1X Assay buffer.</p>	
<p>4. Dilute appropriate conjugate 1/10,000 – 1/100,000 (of a 0.1 mg/mL stock) in freshly prepared conjugate diluent using a minimum of 0.2 mL/cm² of membrane.</p>	<p><i>Example: 2 µL conjugate + 20 mL diluent. Suggested diluents include Detector Block and TBS/PBS-TWEEN. The optimal dilution may vary for different lots of conjugate. Titrate the conjugate to determine the optimal working dilution.</i></p>	<p>10. Gently pipette 0.05 mL/cm² of PhosphaGLO Reserve over the entire membrane. Incubate without rocking for 1 minute.</p>	<p><i>Example: for a 10 x 10 cm blot, use 5 mL of PhosphaGLO Reserve.</i></p> <p><i>The surface tension of the substrate will keep it on the surface of the membrane.</i></p>
<p>5. Incubate blot with diluted conjugate for one hour at room temperature.</p>		<p>11. Lift the membrane with forceps and blot the excess substrate onto a piece of filter paper. Seal the membrane in clear plastic and expose to X-ray film for 10 seconds to 1 minute. Adjust exposure time as needed.</p>	<p><i>Excessive substrate on the blot will contribute to background.</i></p> <p><i>Take caution to ensure the surface of the membrane to which the assay reagents were applied is facing the film. Do not allow the film to get wet, nor move during exposure.</i></p>
<p>6. During the conjugate incubation step, prepare PhosphaGLO Reserve. Prepare 0.05 mL/cm² of membrane to be detected.</p>	<p><i>Allow PhosphaGLO Reserve to come to room temperature prior to use. Cover with foil to minimize light exposure.</i></p>	<p>12. <i>Optional:</i> Chemiluminescent Imager Detection. Incubate the blot for twice the time typically used for film. If the imager provides stacking capabilities, capture exposures at 5 minute intervals for 1 hour to maximize signal. The optimal exposure can be chosen.</p>	<p><i>Optimal exposure time should be determined by the signal to noise ratio and the amount of conjugate used. When using greater amounts of conjugate, 10 seconds may provide acceptable results.</i></p> <p><i>Follow the manufacturer's recommendations regarding the set up and operation of the imager.</i></p>

TROUBLESHOOTING**Problem 1: No Signal**

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Inactive alkaline phosphatase 	Verify enzyme activity by mixing 10 μ L of diluted conjugate with 1 mL of substrate (in a dark room, the substrate should glow).
<ul style="list-style-type: none"> No binding of conjugate to the primary antibody 	Confirm correct specificity of the conjugate for the primary antibody; <i>i.e.</i> no anti-rabbit AP with a mouse primary antibody.
<ul style="list-style-type: none"> No transfer of target to membrane 	Use a protein stain on unblocked membrane to verify attachment of target protein or use a pre-stained protein marker to monitor transfer.
<ul style="list-style-type: none"> Detection of non-blotted side of membrane 	Ensure correct orientation of the membrane during the assay and film exposure.
<ul style="list-style-type: none"> Missed step in procedure 	Review procedure to ensure all steps were followed.

Problem 2: Weak Signal

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Insufficient amount of antibody 	Optimize antibody concentrations. Affinity of the primary antibody may change after proteins are denatured through SDS-PAGE.
<ul style="list-style-type: none"> Insufficient protein loaded or transferred 	Increase the amount of protein loaded onto the gel.
<ul style="list-style-type: none"> Insufficient incubation of primary antibody to target 	Increase the incubation times for weak primary antibodies.
<ul style="list-style-type: none"> Insufficient exposure time 	Increase the time of exposure to film.
<ul style="list-style-type: none"> Excessive washing beyond recommended procedure 	Follow the procedure as written.

Problem 3: Excessive signal, nonspecific bands or general background

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Overexposure of film to signal 	Expose the membrane to film for a shorter period of time.
<ul style="list-style-type: none"> Insufficient blocking or washing 	Increase blocking and washing time or increase number of washes.
<ul style="list-style-type: none"> Excessive antibody used for detection 	Optimize conjugate concentration. Reduce antibody concentrations; optimal conjugate dilution should be 1/10,000 – 1/100,000 of a 0.1 mg/mL stock. OR Decrease the amount of primary antibody.
<ul style="list-style-type: none"> Excessive protein loaded on the gel 	Decrease the amount of protein loaded onto the gel.

Problem 4: Poorly Defined or “Fuzzy” Bands or Dots

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Poor transfer of protein to membrane 	Follow manufacturer’s recommended procedure or contact the manufacturer for additional support regarding the blotting apparatus.
<ul style="list-style-type: none"> Excessive substrate 	Remove excess substrate before exposure of the membrane to film.
<ul style="list-style-type: none"> Ghost images from shifted position of film during development 	Avoid movement of film over membrane during exposure period.
<ul style="list-style-type: none"> Inadequate handling of membranes 	Certain membranes require special handling. Check with the membrane vendor for correct procedures.

Stripping and Reprobing a Western Blot

This protocol is adapted from Kaufmann, *et. al.*¹¹. After performing protein transfer, detection with PhosphaGLO Reserve and film exposure, membranes may be stripped and reprobed with new primary and secondary antibodies.

1. Strip antibodies by incubating blot for 30 - 90 minutes at 70°C in erasure buffer: 2% SDS (w/v), 62.5 mM Tris-HCl (pH 6.8 at 20°C), 100 mM β-mercaptoethanol.
2. Wash 2 times, for 10 minutes each, in TBS: 10 mM Tris-HCl (pH 7.4 at 20°C), 150 mM NaCl.
3. Block for 2.5 hours in Block Solution.
4. Repeat detection procedure.

REFERENCES

1. Kricka, L. J. (1991). Chemiluminescent and Bioluminescent Techniques. *Clin. Chem.* 37(9): 1472 - 1481.
2. Knecht, D.A. and R.L. Dimond (1984). Visualization of Antigenic Proteins on Western Blots. *Anal. Biochem.* 136: 180 - 184.
3. Blake, M.S., et al (1984). A Rapid, Sensitive Method For Detection of Alkaline Phosphatase Conjugated Antibody on Western Blots. *Anal. Biochem.* 136: 175 - 178.
4. Isacsson, V. and G. Wettermark (1974). Chemiluminescence in Analytical Chemistry. *Anal. Chim. Acta.* 68: 339 - 362.
5. Towbin, H., T. Staehelin and J. Gordon (1979). Electrophoretic Transfer of Proteins From Polyacrylamide Gels to Nitrocellulose Sheets: Procedure and Some Applications. *Proc. Natl. Acad. Sci. USA.* 76: 4350 - 4354.
6. Bittner, M., P. Kupferer and C. F. Morris (1980). Electrophoretic Transfer of Proteins and Nucleic Acids From Slab Gels to Diazobenzyloxymethyl Cellulose or Nitrocellulose Sheets. *Anal. Biochem.* 102: 459 - 471.
7. Burnette, W.N. (1980). "Western Blotting": Electrophoretic Transfer of Proteins From Sodium Dodecyl Sulfate-Polyacrylamide Gels to Unmodified Nitrocellulose or Nitrocellulose Sheets. *Anal. Biochem.* 112: 195 - 203.
8. Reinhart, M.P. and D. Malamud (1982). Protein Transfer From Isoelectric Focusing Gels: The Native Blot. *Anal. Biochem.* 123: 229 - 235.
9. Gooderham, K (1983). *Protein Blotting*. In J. Walker and W. Gaastra (eds.), Techniques in Molecular Biology. Croom Helm Ltd. Publishers, London.
10. Southern, E.M. (1975). Detection of Specific Sequences Among DNA Fragments Separated by Gel Electrophoresis. *J. Mol. Biol.* 98: 503 - 517.
11. Kaufmann, Ewing and Shaper (1987). The Erasable Western Blot. *Anal. Biochem.* 161, 89 - 95.

PhosphaGLO and PhosphaGLO Reserve are trademarks of KPL, Inc.

Detector is a trademark of KPL, Inc.

TWEEN is a trademark of ICI Americas, Inc.

BIOMAX is a trademark of Kodak.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However, said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product. No claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims regarding product performance must be made within 30 days following date of delivery.

Limited Use License

The purchase of this product conveys to the buyer the non-transferable right to use the product in research conducted by the buyer. The buyer cannot sell or otherwise transfer this product or materials made by use of this product to a third party or otherwise use this product or materials made with this product for Commercial Purposes without written approval of KPL, Inc. For additional information, please visit www.kpl.com or the KPL catalog. To obtain a license or approval to use this product for purposes other than those permitted above, contact Director of Sales at (301) 948-7755.