

Blocking Solution Concentrate for Immunohistochemistry

Catalog No.
71-00-10

Size
10 mL



DESCRIPTION

Blocking Solution Concentrate for Immunohistochemistry is used to block endogenous peroxidase activity in immunohistochemistry applications.

FORM/STORAGE

Blocking Solution Concentrate is a 10X liquid concentrate. Store tightly capped at 2 - 8°C. Stable for a minimum of one year stored at 2 - 8°C.

USE

Prepare working solution by diluting 1:10 in reagent quality water. Working solution should be prepared fresh before each use.

After tissue sections or other cellular preparations have been fixed or deparaffinized as necessary, immerse specimen in Blocking Solution at working dilution. Suggested incubation times are as follows:

Paraffin Sections:	4 minutes
Frozen Sections, Cytospin Preparations or Blood Films:	45 seconds

After incubation with the blocking solution, rinse specimens in reagent quality water, then soak in pH 7.6 Tris-HCl for 5 minutes before proceeding with the immunoassay.

IMMUNOASSAY KIT REAGENTS

Blocking Solution Concentrate for Immunohistochemistry is a component of the HistoMark kits for Horseradish Peroxidase Substrate Systems:

HistoMark® BLACK	Catalog No. 54-75-00
HistoMark® ORANGE	Catalog No. 54-74-00
HistoMark® TrueBlue	Catalog No. 54-78-00

RELATED PRODUCTS

HistoMark® RED Phosphatase System	Catalog No. 55-69-00
HistoMark® BLUE Phosphatase System	Catalog No. 55-70-00
Universal Block for Immunohistochemistry	Catalog No. 71-00-61

See KPL catalog for a complete list of antibodies, substrates and complete systems for microwell ELISA, membrane blotting and immunohistochemical applications.

PRODUCT SAFETY AND HANDLING

This product is considered non-hazardous as defined by the Hazard Communication Standard (29 CFR 1910.1200). Avoid contact with skin and eyes. In case of contact or spillage, clean with copious amounts of water. Disposal via sanitary sewer.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However, said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product. No claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims regarding product performance must be made within 30 days following date of delivery.

Limited Use License

The purchase of this product conveys to the buyer the non-transferable right to use the product in research conducted by the buyer. The buyer cannot sell or otherwise transfer this product or materials made by use of this product to a third party or otherwise use this product or materials made with this product for Commercial Purposes without written approval of KPL, Inc. For additional information, please visit www.kpl.com or the KPL catalog. To obtain a license or approval to use this product for purposes other than those permitted above, contact Director of Sales at (301) 948-7755.