IVD DATA SHEET

Vimentin

Concentrated Rabbit Monoclonal Antibody

Intended Use:

For in Vitro Diagnostic Use

Epitomics' Rabbit Monoclonal Anti-Human Vimentin, Clone EP21, is intended for use to qualitatively identify Vimentin by light microscopy in sections of formalin-fixed, paraffin-embedded tissue using immunohistochemical detection methodology. Interpretation of any positive or negative staining must be complemented with the evaluation of proper controls and must be made within the context of the patient's clinical history and other diagnostic tests. Evaluation must be performed by a qualified pathologist.

Catalog number	Description	Dilution
AC-0024A*	0.1 ml, concentrated	1:100-1:200
AC-0024B*	0.5 ml, concentrated	1:100-1:200
AC-0024	1 ml, concentrated	1:100-1:200
AC-0024BULK*	2 ml or more, concentrated	1:100-1:200

* Product not available in Canada

Immunogen: A synthetic acetylated peptide corresponding to the C-term of human Vimentin protein was used.

Source: Rabbit Monoclonal Antibody

Clone ID: EP21

Isotype: Rabbit IgG

Application: Immunohistochemistry for formalin-fixed paraffin-embedded tissue

Summary and Explanation:

Vimentin is the most common member of intermediate filament (IF) family and one of the main components in cytoskeleton structure. It is expressed during cell development and differentiation in variety of mesenchymal cells and cell types derived from mesoderm. Vimentin is essential in the role of cell integrity and cytoskeletal stability. The reorganization of Vimentin, similar to all IF proteins, occurs during different stages of the cell cycle and cell signaling by a site-specific phosphorylation (serine and threonine residues). In particular, p21-activated kinase (PAK) phosphorylates at Ser25, Ser38, Ser50, Ser65 and Ser72 which induces Vimentin specific reorganization. During cytokinesis, Vimentin is regulated by Rho-kinase (ROCK) and Aurora B through phosphorylation at Ser38 and Ser72. Also, coordinated by ROCK and Aurora B, Plk1 induces phosphorylation at Ser82 plays an important role in Vimentin segregation. Antibody to Vimentin labels mesenchymal cells. In combination with a panel of antibodies, it is used to identify tumor with mesenchymal origin and malignant melanoma. Additionally, vimentin is a useful control marker for proper tissue processing.

Reagent Provided:

Antibody to Vimentin is affinity purified and diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN₃).

Storage and Stability:

Store at 2-8 °C. Don't use after expiration date provided on the vial. End user must validate any storage conditions other than those specified.

Procedures Recommended:

- 1. Pretreatment:** Epitope retrieval using citrate buffer (catalog #: SP-0001) with a pressure cooker.
- 2. Endogenous peroxidase block:** Block for 10 minutes at room temperature using peroxidase solution (catalog #: SP-0002).
- 3. Protein block:** Block for 10 minutes at room temperature using blocking solution (catalog #: SP-0003).
- 4. Primary antibody:** Incubate for 30 minutes.
- 5. Detection:** Follow instructions from the selected detection system (EpiPrecision™, a Biotin Streptavidin-HRP Detection, catalog #: DK-0001, 0003, or EpiVision™, a Rabbit Polymer Detection, catalog # DK-0002, 0004).

The antibody dilution and protocol may vary depending on the specimen preparation and specific application. Optimal conditions should be determined by the individual laboratory.

Performance Characteristics:

This antibody gives cytoplasm staining in positive cells. The recommended positive controls are Colon for normal tissue and Sarcoma for abnormal tissue.

Limitations:

Immunohistochemistry is a complex process. Variation in tissue selection, tissue processing, antigen retrieval, peroxidase activity, detection systems and improper counterstaining may cause variation in results.

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