IVD DATA SHEET

S100 Beta

Concentrated Rabbit Monoclonal Antibody

Intended Use:

For in Vitro Diagnostic Use

Epitomics' Rabbit Monoclonal Anti-Human S100 Beta, Clone EP32, is intended for use to qualitatively identify S100 Beta by light microscopy in sections of formalin-fixed, paraffin-embedded tissue using immunohistochemical detection methodology. Interpretation of any positive or negative staining must be complemented with the evaluation of proper controls and must be made within the context of the patient's clinical history and other diagnostic tests. Evaluation must be performed by a qualified pathologist.

Catalog number	Description	Dilution
AC-0031A*	0.1 ml, concentrated	1:100-1:200
AC-0031B*	0.5 ml, concentrated	1:100-1:200
AC-0031	1 ml, concentrated	1:100-1:200
AC-0031BULK*	2 ml or more, concentrated	1:100-1:200

* Product not available in Canada

Immunogen: A synthetic peptide corresponding to residues on the C-terminus of human S100 Beta protein

Source: Rabbit Monoclonal Antibody

Clone ID: EP32

Isotype: Rabbit IgG

Application: Immunohistochemistry for formalin-fixed paraffin-embedded tissue

Summary and Explanation:

S100 Beta is a homodimeric member of the S100 superfamily. S100 is a family of Ca²⁺-binding proteins, comprises 19 members that are differentially expressed in a large number of cell types. The protein has been implicated in cellular processes such as cell differentiation and growth. S100 Beta is abundant in glial cells of the central and peripheral nervous system, in melanocytes, chondrocytes, and adipocytes. Antibody to S100 Beta also labels Langerhans cells, histiocytes, epithelial, myoepithelial cells and integrating reticular cells of lymphoid tissue, and tumors originated from these cells. S100 Beta is a useful marker for diagnosis of melanoma, tumors of nerves system.

Reagent Provided:

Antibody to S100 Beta is affinity purified and diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN₃).

Storage and Stability:

Store at 2-8 °C. Don't use after expiration date provided on the vial. End user must validate any storage conditions other than those specified.

Procedures Recommended:

1. Pretreatment: Epitope retrieval using citrate buffer (catalog #: SP-0001) with a pressure cooker.

2. Endogenous peroxidase block: Block for 10 minutes at room temperature using peroxidase solution (catalog #: SP-0002).

3. Protein block: Block for 10 minutes at room temperature using blocking solution (catalog #: SP-0003).

4. Primary antibody: Incubate for 30 minutes.

5. Detection: Follow instructions from the selected detection system (EpiPrecision™, a Biotin Streptavidin-HRP Detection, catalog #: DK-0001, 0003, or EpiVision™, a Rabbit Polymer Detection, catalog # DK-0002, 0004).

The antibody dilution and protocol may vary depending on the specimen preparation and specific application. Optimal conditions should be determined by the individual laboratory.

Performance Characteristics:

This antibody gives cytoplasmic staining in positive cells. The recommended positive controls are skin for normal tissue and melanoma for abnormal tissue.

Limitations:

Immunohistochemistry is a complex process. Variation in tissue selection, tissue processing, antigen retrieval, peroxidase activity, detection systems and improper counterstaining may cause variation in results.

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