IVD DATA SHEET

Androgen Receptor

Concentrated Rabbit Monoclonal Antibody

Intended Use:

For in Vitro Diagnostic Use

Epitomics' Rabbit Monoclonal Anti-Human Androgen receptor (AR), Clone EP120, is intended for use to qualitatively identify Androgen Receptor by light microscopy in sections of formalin-fixed, paraffin-embedded tissue using immunohistochemical detection methodology. Interpretation of any positive or negative staining must be complemented with the evaluation of proper controls and must be made within the context of the patient's clinical history and other diagnostic tests. Evaluation must be performed by a qualified pathologist.

Catalog number	Description	Dilution
AC-0071A	0.1 ml, concentrated	1:100-1:200
AC-0071B	0.5 ml, concentrated	1:100-1:200
AC-0071	1 ml, concentrated	1:100-1:200
AC-0071BULK	2 ml or more, concentrated	1:100-1:200

Immunogen: A synthetic peptide corresponding to residues in the N-terminus of human Androgen Receptor protein

Source: Rabbit Monoclonal Antibody

Clone ID: EP120

Isotype: Rabbit IgG

Application: Immunohistochemistry for formalin-fixed paraffin-embedded tissue

Summary and Explanation:

Androgen receptor (AR) is a member of the steroid receptor superfamily that is essential for the growth of prostate cancer cells. It has been reported that tyrosine phosphorylation of AR is induced by growth factors and elevated in hormone-refractory prostate tumors. Data suggest that growth factors and their downstream tyrosine kinases, which are elevated during hormone-ablation therapy, can induce tyrosine phosphorylation of AR. Such modification may be important for prostate tumor growth under androgen-depleted conditions. Cellular signaling occurs following androgen binding to the AR and translocation to the nucleus. This activated complex associates with androgen-responsive elements contained in the DNA sequence of target genes, affecting the transcriptional activity of these genes.

AR antibody labels epithelial cells and stromal cells in normal prostate. AR reactivity is also found in other types of cells, including epithelial cells of the breast and hepatocytes. In prostate cancer, AR expression is maintained throughout cancer progression. Immunohistochemistry of AR is useful for the evaluation of prostate cancer AR in routinely processed tissues. The majority of androgen-independent hormone refractory prostate cancers express AR. Transcriptional activation of AR is involved in refractory anti-androgen therapy. Androgen receptor expression has been helpful in predicting the clinical response to anti-androgenic treatment.

Reagent Provided:

Antibody to Androgen Receptor is affinity purified and diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN₃).

Storage and Stability:

Store at 2-8 °C. Don't use after expiration date provided on the vial. End user must validate any storage conditions other than those specified.

Procedures Recommended:

1. Pretreatment: Epitope retrieval using citrate buffer (catalog #: SP-0001) with a pressure cooker.

2. Endogenous peroxidase block: Block for 10 minutes at room temperature using peroxidase solution (catalog #: SP-0002).

3. Protein block: Block for 10 minutes at room temperature using blocking solution (catalog #: SP-0003).

4. Primary antibody: Incubate for 30 minutes.

5. Detection: Follow instructions from the selected detection system (EpiPrecision™, a Biotin Streptavidin-HRP Detection, catalog #: DK-0001, 0003, or EpiVision™, a Rabbit Polymer Detection, catalog # DK-0002, 0004).

The antibody dilution and protocol may vary depending on the specimen preparation and specific application. Optimal conditions should be determined by the individual laboratory.

Performance Characteristics:

This antibody gives nuclear staining in positive cells. The recommended positive controls are Prostate for normal tissue and Prostate cancer for abnormal tissue.

Limitations:

Immunohistochemistry is a complex process. Variation in tissue selection, tissue processing, antigen retrieval, peroxidase activity, detection systems and improper counterstaining may cause variation in results.

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