IVD DATA SHEET

Topoisomerase II alpha

Concentrated Rabbit Monoclonal Antibody

Intended Use:

For in Vitro Diagnostic Use (For distribution outside United States only)

Epitomics' Rabbit Monoclonal Anti-Human Topoisomerase II alpha (TOP2A), Clone EP93, is intended for use to qualitatively identify Topoisomerase II alpha by light microscopy in sections of formalin-fixed, paraffin-embedded tissue using immunohistochemical detection methodology. Interpretation of any positive or negative staining must be complemented with the evaluation of proper controls and must be made within the context of the patient's clinical history and other diagnostic tests. Evaluation must be performed by a qualified pathologist.

Catalog number	Description	Dilution
AC-0089EUA	0.1 ml, concentrated	1:100-1:200
AC-0089EUB	0.5 ml, concentrated	1:100-1:200
AC-0089EU	1 ml, concentrated	1:100-1:200
AC-0089EUBULK	2 ml or more, concentrated	1:100-1:200

Immunogen: A synthetic peptide corresponding to C-terminal residues of human Topoisomerase II alpha (TOP2A) protein

Source: Rabbit Monoclonal Antibody

Clone ID: EP93

Isotype: Rabbit IgG

Application: Immunohistochemistry for formalin-fixed paraffin-embedded tissue

Summary and Explanation:

DNA topoisomerase II alpha (Topo IIα) is a 170 kDa ubiquitous nuclear enzyme belonging to the topo family, which regulates the topological states of DNA. Topo II is required in chromatin condensation and segregation during mitosis. Topo II α is cell cycle regulated and its level peaks between G2 to M phase. It has been linked to cell proliferation and it may be the main isoform of Topo II involved mitotic processes. Topo II α passes one strand of DNA through a reversible break in a second DNA strand, which catalyzes the topological isomerization of DNA during cell cycle. Topo II α overexpression has been linked to a number of human malignancies and is the target for many chemotherapeutic agents. The majority of anticancer drugs targeting Topo IIα initiate apoptosis by stabilizing the covalent complex formed between DNA and Topo IIα.

Reagent Provided:

Antibody to TOP2A is affinity purified and diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN₃).

Storage and Stability:

Store at 2-8 °C. Don't use after expiration date provided on the vial. End user must validate any storage conditions other than those specified.

Procedures Recommended:

- 1. Pretreatment:** Epitope retrieval using citrate buffer (catalog #: SP-0001) with a pressure cooker.
- 2. Endogenous peroxidase block:** Block for 10 minutes at room temperature using peroxidase solution (catalog #: SP-0002).
- 3. Protein block:** Block for 10 minutes at room temperature using blocking solution (catalog #: SP-0003).
- 4. Primary antibody:** Incubate for 30 minutes.
- 5. Detection:** Follow instructions from the selected detection system (EpiPrecision™, a Biotin Streptavidin-HRP Detection, catalog #: DK-0001, 0003, or EpiVision™, a Rabbit Polymer Detection, catalog # DK-0002, 0004).

The antibody dilution and protocol may vary depending on the specimen preparation and specific application. Optimal conditions should be determined by the individual laboratory.

Performance Characteristics:

This antibody gives nuclear/cytoplasmic staining in positive cells. The recommended positive controls are Tonsil for normal tissue and Gastric carcinoma for abnormal tissue.

Limitations:

Immunohistochemistry is a complex process. Variation in tissue selection, tissue processing, antigen retrieval, peroxidase activity, detection systems and improper counterstaining may cause variation in results.

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