

UnoVue™ Mouse AP Detection Reagent

Catalog No. MU-AP100, MU-AP1000
Document #: DS-6029-A
Effective Date: 05/01/2015

Intended Use For In Vitro Diagnostic Use

UnoVue Mouse AP Detection Reagent is a non-biotin one-step detection reagent suitable for demonstrating antigens in formalin-fixed paraffin-embedded tissues and frozen sections. The UnoVue detection reagents may also be used with blood smears, cytospins, and cell preparations.

UnoVue detection reagents have been developed by directly labeling anti-mouse immunoglobulins with enzymes using a proprietary tandem hyperlabelling technology. This ensures consistent and reproducible immunodetection of mouse antibodies against nuclear, cytoplasmic and membrane antigens in different types of tissues. The single step UnoVue Detection Reagent enables faster staining procedures than traditional two-step methods using biotin and avidin/streptavidin conjugates, with significantly lower background.

UnoVue Mouse AP Detection Reagent is suitable for use with mouse all mouse antibodies, both monoclonal and polyclonal. The reagents can be used for manual staining or with automated staining instruments and are well suited for multiplex immunohistochemical staining assays.

Kit Contents

Description	Catalog #	Volume
UnoVue Mouse AP Detection Reagent	MU-AP100	10 ml (100 Tests)
	MU-AP1000	100 ml (1000 Tests)

Storage and Handling Store at 2°-8°C away from light. Do not use product after the expiration date printed on vial. If reagents are stored under conditions other than those specified here, they must be verified by the user. Diluted reagents should be used promptly.

Stability 12-18 months (see expiration date on reagent bottles)

Composition UnoVue Mouse AP ready to use reagent is formulated without azide or thimerosal preservatives.

Material Required But Not Provided Some of the reagents and materials required for IHC are not provided. Pretreatment reagents, detection systems, control reagents and other ancillary reagents are available from Diagnostic BioSystems. Please refer to the Diagnostic BioSystems website at www.dbiosys.com.

Precautions

- Wear appropriate personal protective apparel. Avoid contact with clothes and exposed skin. In case of accidental skin exposure, flush with water immediately. Consult a physician if required.
- Interpretation of the results is the sole responsibility of the user.

Troubleshooting If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem is suspected, contact Diagnostic BioSystems Technical Support at (925) 484-3350, extension 2 or techsupport@dbiosys.com.

Recommended Staining Protocol

- Paraffin embedded tissue sections must be deparaffinized with xylene or dewaxing agent and rehydrated with a graded series of ethanol and water washes before staining. Follow the standard dewaxing and rehydration protocol used in your lab.
- The investigator needs to optimize the dilution and incubation times for primary antibodies.
- Each immunostaining run should include known positive and negative controls to assure proper functioning of the staining system and aid in valid interpretation of the results.

Typical controls:



Diagnostic BioSystems
 6616 Owens Drive
 Pleasanton, CA, 94588
 Tel: (925) 484 3350
www.dbiosys.com



Emergo Europe
 Molenstraat 15
 2513 BH, The Hague
 The Netherlands
 Tel: (31) (0) 70 345
 8570

Positive Control: A tissue known to contain the desired antigen, which has yielded positive staining in the past.

Negative Controls:

Reagent Controls

- A. Substitute normal non-immune serum from the same host animal as the primary antibody (e.g. if using mouse monoclonal primary antibodies, use mouse non-immune serum).
- B. Substitute matching host species isotype control for primary antibody
- C. Use antigen-adsorbed primary antibody (i.e. antibody reagent which has been adsorbed with the target antigen to remove specific antibody)

Tissue control – A tissue known to *not* contain the desired antigen.

- 4. Consult the primary antibody supplier for recommended for antigen recovery treatments. Perform epitope recovery pretreatments before starting the staining procedure.
- 5. Once the slide treatment has been started, **DO NOT** let tissues or specimens dry. This can cause undesirable background or artifacts.
- 6. After the Primary Antibody application wash slides with Immuno Wash buffer (DBS cat# K 005).
- 7. Apply UnoVue Mouse AP reagent (Make sure tissue sections/smears completely covered with the reagent) and incubate for 20 min at room temperature.
- 8. Wash slides with Immuno Wash buffer then proceed to chromogen labeling and counterstain steps.



Diagnostic BioSystems
6616 Owens Drive
Pleasanton, CA, 94588
Tel: (925) 484 3350
www.dbiosys.com



Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: (31) (0) 70 345
8570