

medac



Quality in Control

# PD-L1 Analyte Control<sup>DR</sup>



Theaterstraße 6  
D-22880 Wedel

Telefon 04103/ 8006-342  
Telefax 04103/ 8006-359


[www.medac-diagnostika.de](http://www.medac-diagnostika.de)

Product Codes: HCL019, HCL020 and HCL021

# Contents

PD-L1 Analyte Control <sup>DR</sup>	2
What is PD-L1?	3
The Role of PD-L1 in Cancer	3
PD-L1 Assessment	4
PD-L1 Product Details	5
PD-L1 Analyte Control Staining	6


HistoCyte Laboratories Ltd is based in the heart of the Newcastle University campus. Started in 2014 by scientists with a combined experience of over 30 years in the development of reagents for immunohistochemistry and in-situ hybridization. Collaborating with pathologists locally and globally, HistoCyte Laboratories Ltd is developing a range of cost effective products designed to help scientists to maintain and develop the quality of assays within their laboratory.




# PD-L1 Analyte Control<sup>DR</sup>

**PD-L1 Analyte Control<sup>DR</sup>** is part of the *Dynamic Range* of HistoCyte Products. When a Dynamic Range of controls demonstrating the sensitivity of an assay is required, the **PD-L1 Analyte Control<sup>DR</sup>** is ideal. This product contains four cells of varying expression including a negative control cell.

**PD-L1 Analyte Control<sup>DR</sup>** is available as pre-cut slides (2 slide and 5 slides) and cell microarray blocks.



Format	Product Code
2 Slide	HCL019
5 Slide	HCL020
Block	HCL021



# What is PD-L1?

Programmed death ligand 1 (PD-L1) is a 40kD type 1 transmembrane protein. Synonyms include:

- CD274
- B7 homolog 1 (B7-H1)

PD-L1 is a checkpoint regulator in immune cells<sup>1</sup>, it is expressed on immune or non-hematopoietic cells<sup>2</sup>. Expression of the protein is seen during pregnancy where it has a role in suppressing the immune system. PD-L1 induces an inhibitory signal in activated T-cells and promotes T-cell apoptosis<sup>2</sup>.

## The Role of PD-L1 in Cancer

PD-L1 has been observed to be over expressed in a number of different cancer types and is believed to be a potential means by which the cancer cells can evade the immune system. Overexpression of PD-L1 correlates with poor disease outcomes<sup>3</sup>. The expression of PD-L1 within cancer is not restricted to a single type of cancer and as such it has become a target for anti-cancer drug development. Currently there are a number of anti-PD-L1 clinical trials ongoing, focusing on the following tumour types:

- Lung cancer
- Bladder cancer
- Kidney cancer
- Haematological cancer
- Breast cancer
- Colorectal cancer
- Melanoma
- Solid tumours

1. Dong H et al Nat Med 1999 5 1365-1369

2. Shi L et al J Hematol Oncol. 2013; 6: 74.

3. Ohaegbulam K, et al Cell 2015 21, Issue 1, p24–33

The diagram below (Figure 1.) illustrates the interaction between the tumour cells and the immune system, whereby the anti-PD-L1 antibody blocks the ability of the ligand to bind with the PD-1 receptor. Thus preventing the inhibitory feedback that would otherwise be stimulated.




Figure 1. The human leukocyte antigen (HLA) on the tumour cell presents tumour protein which is detected through the T cell receptor (TCR). Upon recognising “tumour protein” the T cell initiates a cytotoxicity event, which would otherwise be inhibited by the interaction between PD-L1 and PD-1 on the T cell.

## PD-L1 Assessment

A number of different methods are used to measure PD-L1 expression, these include molecular methods such as:

- Real-time polymerase chain reaction (PCR) using products such as TaqMan<sup>®</sup> gene expression assay from ThermoFisher.

- Fluorescence in situ hybridisation (FISH) probes for the detection of PD-L1 DNA.
- Advanced Cell Diagnostics provide an RNAscope product for the detection of PD-L1 mRNA.

A number of antibodies are available for the immunohistochemical detection of PD-L1, these include clones:

- E1L3N (Cell Signalling Technology)
- SP263 (Ventana, Roche)
- SP142 (Spring Bioscience)
- 28-8 (Dako, Agilent).
- 22C3 (Dako, Agilent).

## PD-L1 Product Details


The product consists of four different cell lines with PD-L1 expression levels of high, medium, low and negative. The product was developed using the E1L3N PD-L1 antibody and SP263. However, it has also been independently tested at different laboratories using different PD-L1 antibodies with different protocols. In all cases the HistoCyte PD-L1 control provided the same high, medium, low and zero expression range regardless of assay employed. All HistoCyte products are designed to be suitable for FISH testing. All the cells are amplified for PD-L1 except the negative cell line.

# PD-L1 Analyte Control<sup>DR</sup> staining

## Cell Signalling Technologies Clone: SP263


Performed on the Roche Ventana Benchmark Ultra™

T cell non-Hodgkin  
Lymphoma




**High expression:** Strong staining in majority of cells.

Fibrosarcoma




**Medium expression:** Convincing staining in majority of cells. Some strong staining.

Osteosarcoma



**Low expression:** Faint staining in majority of cells. Occasional strong staining.

Breast Ductal  
Carcinoma



**Negative expression:** Absence of any genuine staining.

**Additional assessment:** mRNA levels by QT-PCR correlated to the IHC results.



# Also Available from HistoCyte Laboratories Ltd

Product Name	Format	Code
<b>HPV/p16 Analyte Control<sup>DR</sup></b> (Four core with dynamic range of HPV gene copies)	Slide(2)	HCL001
	Slide(5)	HCL002
	Block	HCL003
<b>HPV/p16 Analyte Control</b> (Three core with standard range of HPV gene copies)	Slide(2)	HCL004
	Slide(5)	HCL005
	Block	HCL006
<b>ALK-Lung Analyte Control</b> (Two core positive and negative for the EML4-ALK translocation)	Slide(2)	HCL007
	Slide(5)	HCL008
	Block	HCL009
<b>ALK-Lymphoma Analyte Control</b> (Two core positive and negative for the NPM-ALK translocation)	Slide(2)	HCL010
	Slide(5)	HCL011
	Block	HCL012
<b>Breast Analyte Control</b> (Two cores, one positive for Her2, ER and PR. The other negative)	Slide(2)	HCL013
	Slide(5)	HCL014
	Block	HCL015
<b>Breast Analyte Control<sup>DR</sup></b> (Five cores with a dynamic range of expression of Her2, ER and PR. Including negative control)	Slide(2)	HCL016
	Slide(5)	HCL017
	Block	HCL018
<b>PD-L1 Analyte Control<sup>DR</sup></b> (4 core with a dynamic range of expression of PD-L1)	Slide(2)	HCL019
	Slide(5)	HCL020
	Block	HCL021
<b>ROS1 Analyte Control</b> (Two cores positive and negative for ROS1 translocation)	Slide(2)	HCL022
	Slide(5)	HCL023
	Block	HCL024
<b>Sienna Cancer Diagnostics hTERT assay.</b> 1ml of anti-hTERT mouse mAb. <i>(Available UK &amp; Ireland Only)</i>	1ml	HCL025

Your local distributor:



Theaterstraße 6      Telefon 04103/ 8006-342  
D-22880 Wedel      Telefax 04103/ 8006-359  
[www.medac-diagnostika.de](http://www.medac-diagnostika.de)



Quality in Control