

Changes in the Instruction for Use of new CD-ROM 05/2014

All medac Kits:

NEW	OLD
2.1. Distilled or deionised water.	2.1. Water for injection (H ₂ O redist.). Use of deionised water can disturb the test procedure.
4.1. The test is suitable for serum samples and EDTA-plasma. (excluded HAMA)	4.1. The test is suitable for serum samples.

Reference numbers:

430-TMB/ 431-TMB /432-TMB
480-TMB/ 485-TMB /490-TMB
497-TMB/ 498-TMB

NEW	OLD
KIT CONTENTS 6. CON Conjugate: 3 vials with 4.5 ml each, goat anti-human IgX	KIT CONTENTS 6. CON Conjugate: 4 vials with 4.5 ml each, goat anti-human IgX

Reference number:

103-PKS

NEW	OLD
KIT CONTENTS 10. BLOCK removed	KIT CONTENTS 10. BLOCK Blocking reagent: 1 vial with 10 ml