


High performance **in harmony**



Product Specification
BioCLIA® 500
www.bioclia.de

BioCLIA® 500

Fully automated Chemiluminescence Analyzer for Autoimmune and Allergy Diagnostics[#]

System Specifications:

- 72 test results / hour
- 48 autoimmune parameters CE-certified (see product list)
- 130 cuvettes load capacity, refill anytime
- Built-in PC touch screen
- 11 autoimmune or 33 allergy onboard reagent positions
- Temperature control 2 - 8 °C
- 24 sample positions incl. 1 STAT position, random-access
- Continued loading of samples and reagents
- Barcoded QC and calibrators
- Stability of opened reagents and calibrators: 28 days
- Bi-directional HL7- LIS connection
- Dimension (L×W×H): 500mm × 645mm × 648mm
- Weight: 89kg
- Power Supply: 100 - 240V~, 50Hz / 60Hz ± 1Hz, 800VA

[#]30 allergenes CE-certified
available autumn 2022

Diagnostics in a new dimension -
faster, quieter, more efficient and compact!

BioCLIA® 500 – Chemiluminescence Analyzer

