

BioCLIA® 500

Your **compact** solution for the
Autoimmune and Allergy diagnostics[#]

- ✓ Innovative, fully automated CLIA technology
- ✓ Unsurpassed speed, high throughput
- ✓ Accurate results through better discrimination
- ✓ Comprehensive parameter portfolio
- ✓ Easy, intuitive handling
- ✓ Random access, continuous load



Specific Controls
and Calibrators for
each kit available.



Indication list
BioCLIA® 500
www.bioclia.de

BioCLIA® 500

Fully Automated Chemiluminescence Analyzer

Available parameter portfolio for the Autoimmune Disease Areas:**ENA**

- CENP-B
- DFS-70*
- dsDNA
- Histone
- Jo-1
- nRNP/Sm
- Nucleosome
- PCNA
- PM-Scl
- Ribosomal P
- Ro52
- Ro60
- Scl-70
- Sm
- SSB (La)
- ENA-Screen¹
- CTD-Screen²

Liver

- AMA-M2
- F-actin*
- Gp210
- LC-1
- LKM-1
- SLA/LP
- sp100

Celiac disease

- DGP IgA
- DGP IgG
- tTG IgA
- tTG IgG

Antiphospholipid syndrome

- Cardiolipin IgA
- Cardiolipin IgG
- Cardiolipin IgM
- Cardiolipin Screen³
- β2 Glycoprotein 1 IgA
- β2 Glycoprotein 1 IgG
- β2 Glycoprotein 1 IgM
- β2 Glycoprotein 1 Screen⁴

Vasculitis

- GBM
- MPO
- PR3

Rheumatoid arthritis

- CCP
- RF IgA
- RF IgG
- RF IgM
- RF Screen⁵

Diabetes

- IA2
- IAA
- ICA

#30 allergenes CE-certified
available autumn 2022

**All products are CE-certified -
further parameters in preparation**

¹ CENP-B, Jo-1, nRNP/Sm, Ro52, Ro60, Scl-70, Sm, SSB (La);

² AMA-M2, dsDNA, Histone, nRNP/Sm, Nucleosome, Rib-P, PCNA, PM-Scl, Ro52, Ro60, Scl-70, Sm, SSB (La), Jo-1, CENP-B;

³ Cardiolipin IgA/IgG/IgM;

⁴ β2 Glycoprotein 1 IgA/IgG/IgM;

⁵ RF IgA/IgG/IgM

* in preparation



Diagnostics in a new dimension -
faster, quieter, more efficient and compact!

BioCLIA® 500 – Chemiluminescence Analyzer