

**MASTER DIAGNÓSTICA, S.L.**

Avda. Conocimiento, 100

P.T. Ciencias de la Salud

18016 GRANADA


www.masterdiagnostica.com

Fax: 958 27.14.34

Tlf: 958 27.14.49

ANTIBODY DILUENT

PRESENTATION

Antibody diluent, ready-to-use . The presentations for this product are as follow:

Referencies	Volume (ml)
MAD-004048R	50
MAD-004048R-100	100
MAD-004048R-250	250

CONTENTS, DESCRIPTION AND APLICATION:

This product contains 1% Bovine Serum Albumin in phosphate buffer, pH 7.4 and 0.05% sodium Azide.

The Antibody Diluent is prepared for use in immunohistochemical procedures. Monoclonal and polyclonal antibodies can be diluted directly in this reagent and can be stored for up to 18 months when refrigerated (extremely high dilutions will result in reduced shelf life; shelf time should be determined by the individual laboratory).

This product does not contain enzymes, steroids, haptens, antigens or other substances that might interfere or potentially change its stability or function.

Note: Please use protection material during it use

SAFETY RECOMMENDATIONS

This product is intended for laboratory professional use only. The product is NOT intended to be used as a drug or for domestic purposes. The current version of the Safety Data Sheet for this product can be downloaded by searching the reference number at www.vitro.bio or can be requested at regulatory.md@vitro.bio.

STORAGE:

Store at 2-8°C. Please do not freeze