

## Information about Compilation of Safety Data Sheets

**Product:** all Antibodies and Antibody conjugates  
**Hersteller:** SeraCare Life Sciences  
910 Clopper Road, Gaithersburg, Maryland 20878 USA  
**Article numbers:** see attachment

All Antibodies and antibody conjugates listed in the attached table had been classified according to Regulation (EC) No. 1907/2006 (REACH) and to Regulation (EC) No. 1272/2008 (GHS/CLP).

None of the Antibody and the antibody conjugates had been classified as hazardous products as defined by the aforementioned Regulations. Therefore, no Safety Data Sheets had been compiled.

When handling the Antibodies and Antibody conjugates, the general and specific safety precautions for laboratories and the corresponding directives have to be followed.

Date of issue: February 2020

**medac**

Product Safety  
Safety, Health & Environment  
Phone: +49 (0)4103 8006 0  
Fax: +49 (0)4103 8006 100  
E-mail: productsafety@medac.de

Theaterstrasse 6  
22880 Wedel  
Germany

Geschäftsführer:  
Jens Denker  
Jörg Hans  
Nikolaus Graf Stolberg  
Heiner Will

Registergericht:  
Pinneberg HRB 12042 PI

[Attachment](#)

| Article Number | Description  |
|----------------|--|
| 5210-0154      | Goat Anti-Human IgA (alpha)  |
| 5210-0155      | Goat Anti-Human IgG (gamma)  |
| 5210-0157      | Goat Anti-Human IgM (my)   |
| 5210-0158      | Goat Anti-Human IgE (epsilon)                                      |
| 5210-0159      | Goat Anti-Human IgG (H+L)  |
| 5210-0160      | Goat Anti-Human IgA+IgG+IgM (H+L)                                  |
| 5210-0164      | Goat Anti-Human IgA + IgG + IgM (H+L), Mouse Serum Adsorbed        |
| 5210-0165      | Goat Anti-Human IgG (Fc)   |
| 5210-0168      | Goat Anti-Bovine IgM (my)  |
| 5210-0169      | Goat Anti-Bovine IgG (H+L)   |
| 5210-0170      | Rabbit Anti-Goat IgG (H+L)   |
| 5210-0171      | Goat Anti-Swine IgG (gamma)  |
| 5210-0172      | Goat Anti-Swine IgM (my)   |
| 5210-0173      | Goat Anti-Swine IgG (H+L)  |
| 5210-0174      | Goat Anti-Rabbit IgG (H+L)   |
| 5210-0176      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                   |
| 5210-0177      | Goat Anti-Rat IgG (gamma)  |
| 5210-0178      | Goat Anti-Rat IgM (my)   |
| 5210-0179      | Goat Anti-Rat IgG (H+L)  |
| 5210-0180      | Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed                      |
| 5210-0181      | Goat Anti-Guinea Pig IgG (H+L)                                     |
| 5210-0182      | Goat Anti-Mouse IgA (alpha)  |
| 5210-0183      | Goat Anti-Mouse IgG (gamma)  |
| 5210-0184      | Goat Anti-Mouse IgM (my), Human Serum Adsorbed                     |
| 5210-0185      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5210-0187      | Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed        |
| 5210-0188      | Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed              |
| 5210-0191      | Goat Anti-Mouse IgG (H+L), Rat Serum Adsorbed                      |
| 5210-0192      | Goat Anti-Mouse IgG (H+L), Rabbit Serum Adsorbed                   |
| 5210-0193      | Goat Anti-Dog IgG (gamma)  |
| 5210-0194      | Goat Anti-Dog IgM (my)   |
| 5210-0195      | Goat Anti-Dog IgG (H+L)  |
| 5210-0197      | Goat Anti-Cat IgM (my)   |
| 5210-0198      | Goat Anti-Cat IgG (H+L)  |
| 5210-0199      | Goat Anti-Horse IgG (gamma)  |
| 5210-0200      | Goat Anti-Horse IgM (my)   |
| 5210-0201      | Goat Anti-Horse IgG (H+L)  |
| 5210-0203      | Rabbit Anti-Sheep IgG (gamma)                                      |
| 5210-0204      | Rabbit Anti-Sheep IgG (H+L)  |
| 5210-0205      | Goat Anti-Chicken IgG (H+L)  |
| 5210-0206      | Goat Anti-Duck IgG (H+L)   |
| 5210-0207      | Goat Anti-Turkey IgG (H+L)   |
| 5210-0209      | Rabbit Anti-Deer IgG (H+L)   |
| 5210-0210      | Goat Anti-Peromyscus leucopus IgG (H+L)                            |
| 5210-0211      | Goat Anti-FITC   |
| 5210-0212      | Goat Anti-PE   |
| 5210-0213      | Goat Anti-APC  |
| 5210-0214      | Goat Anti-PerCP  |
| 5210-0215      | Goat Anti-Monkey IgA (alpha), liquid                               |
| 5210-0216      | Goat Anti-Monkey IgG (gamma), liquid                               |
| 5210-0217      | Goat Anti-Monkey IgM (my), liquid                                  |
| 5210-0297      | Goat Anti-Human IgG (gamma), Fab2 Fragment                         |
| 5210-0298      | Goat Anti-Human IgM (my), F(ab)2 Fragment                          |
| 5210-0299      | Goat Anti-Human IgG (H+L), F(ab)2 Fragment                         |
| 5210-0300      | Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed  |
| 5210-0302      | Goat Anti-Mouse IgG (gamma), F(ab)2 Fragment, Human Serum Adsorbed |
| 5210-0303      | Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed   |
| 5220-0277      | Goat Anti-Human IgG (H+L)  |
| 5220-0278      | Goat Anti-Human IgA+IgG+IgM (H+L), Mouse Serum Adsorbed            |
| 5220-0279      | Goat Anti-Human IgG (Fc)   |
| 5220-0280      | Goat Anti-Bovine IgM (my)  |
| 5220-0281      | Goat Anti-Swine IgG (gamma)  |
| 5220-0282      | Goat Anti-Swine IgM (my)   |
| 5220-0283      | Goat Anti-Rabbit IgG (H+L)   |
| 5220-0284      | Goat Anti-Rat IgG (gamma)  |
| 5220-0285      | Goat Anti-Rat IgM (my)   |
| 5220-0286      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5220-0287      | Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed              |
| 5220-0288      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5220-0289      | Goat Anti-Dog IgG (gamma)  |
| 5220-0290      | Goat Anti-Dog IgM (my)   |
| 5220-0291      | Goat Anti-Cat IgG (gamma)  |
| 5220-0292      | Goat Anti-Cat IgM (my)   |

| Article Number | Description   |
|----------------|---|
| 5220-0293      | Goat Anti-Horse IgG (gamma)                                 |
| 5220-0294      | Goat Anti-Horse IgM (my)                                    |
| 5220-0295      | Rabbit Anti-Sheep IgG (gamma)                               |
| 5220-0296      | Goat Anti-Duck IgG (H+L)                                    |
| 5220-0298      | Rabbit Anti-Deer IgG (H+L)                                  |
| 5220-0299      | Goat Anti-FITC  |
| 5220-0300      | Goat Anti-PE  |
| 5220-0301      | Goat Anti-APC   |
| 5220-0302      | Goat Anti-PerCP   |
| 5220-0303      | Goat Anti-Human IgG (H+L)                                   |
| 5220-0304      | Goat Anti-Human IgG (Fc)                                    |
| 5220-0305      | Goat Anti-Bovine IgM (my)                                   |
| 5220-0306      | Goat Anti-Swine IgG (gamma)                                 |
| 5220-0307      | Goat Anti-Swine IgM (my)                                    |
| 5220-0308      | Goat Anti-Rabbit IgG (H+L)                                  |
| 5220-0309      | Goat Anti-Rat IgG (gamma)                                   |
| 5220-0310      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed             |
| 5220-0311      | Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed       |
| 5220-0312      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed    |
| 5220-0314      | Goat Anti-Dog IgM (my)                                      |
| 5220-0315      | Goat Anti-Cat IgG (gamma)                                   |
| 5220-0316      | Goat Anti-Horse IgG (gamma)                                 |
| 5220-0317      | Goat Anti-Horse IgM (my)                                    |
| 5220-0318      | Rabbit Anti-Sheep IgG (gamma)                               |
| 5220-0320      | Goat Anti-Ferret IgG (H+L)                                  |
| 5220-0322      | Goat Anti-Dog IgG (gamma)                                   |
| 5220-0323      | Goat Anti-FITC  |
| 5220-0324      | Goat Anti-PE  |
| 5220-0325      | Goat Anti-APC   |
| 5220-0326      | Goat Anti-PerCP   |
| 5220-0327      | Goat Anti-Human IgG (gamma)                                 |
| 5220-0328      | Goat Anti-Human IgM (my)                                    |
| 5220-0329      | Goat Anti-Human IgE (epsilon)                               |
| 5220-0330      | Goat Anti-Human IgG (H+L)                                   |
| 5220-0331      | Goat Anti-Human IgA + IgG + IgM (H+L)                       |
| 5220-0332      | Goat Anti-Monkey IgA (alpha)                                |
| 5220-0333      | Goat Anti-Monkey IgG (gamma)                                |
| 5220-0334      | Goat Anti-Monkey IgM (my)                                   |
| 5220-0335      | Goat Anti-Rabbit IgG (H+L), Extra Serum Adsorbed            |
| 5220-0336      | Goat Anti-Rabbit IgG (H+L)                                  |
| 5220-0337      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed            |
| 5220-0338      | Goat Anti-Mouse IgG (H+L), Extra Serum Adsorbed             |
| 5220-0339      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed           |
| 5220-0340      | Goat Anti-Mouse IgM (my), Human Serum Adsorbed              |
| 5220-0341      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed             |
| 5220-0342      | Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed |
| 5220-0343      | Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed       |
| 5220-0347      | Goat Anti-Human IgA (alpha)                                 |
| 5220-0348      | Goat Anti-Human IgG (gamma)                                 |
| 5220-0349      | Goat Anti-Human IgM (my)                                    |
| 5220-0350      | Goat Anti-Human IgE (epsilon)                               |
| 5220-0351      | Goat Anti-Human IgG (H+L)                                   |
| 5220-0352      | Goat Anti-Human IgA + IgG + IgM (H+L)                       |
| 5220-0353      | Goat Anti-Rabbit IgG (H+L)                                  |
| 5220-0354      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed            |
| 5220-0355      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed           |
| 5220-0356      | Goat Anti-Mouse IgM (my), Human Serum Adsorbed              |
| 5220-0357      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed             |
| 5220-0358      | Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed |
| 5220-0359      | Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed       |
| 5220-0360      | Goat Anti-Human IgA (alpha)                                 |
| 5220-0361      | Goat Anti-Bovine IgG (H+L)                                  |
| 5220-0362      | Rabbit Anti-Goat IgG (H+L)                                  |
| 5220-0363      | Goat Anti-Swine IgG (H+L)                                   |
| 5220-0364      | Goat Anti-Rat IgG (H+L)                                     |
| 5220-0365      | Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed               |
| 5220-0366      | Goat Anti-Guinea Pig IgG (H+L)                              |
| 5220-0367      | Goat Anti-Mouse IgA (alpha), Human Serum Adsorbed           |
| 5220-0368      | Goat Anti-Dog IgG (H+L)                                     |
| 5220-0369      | Goat Anti-Cat IgG (H+L)                                     |
| 5220-0370      | Goat Anti-Horse IgG (H+L)                                   |
| 5220-0371      | Goat Anti-Hamster IgG (H+L)                                 |
| 5220-0372      | Rabbit Anti-Sheep IgG (H+L)                                 |

| Article Number | Description  |
|----------------|--|
| 5220-0373      | Goat Anti-Chicken IgG (H+L)  |
| 5220-0374      | Goat Anti-Turkey IgG (H+L)   |
| 5220-0375      | Goat Anti-Peromyscus leucopus IgG (H+L)                            |
| 5220-0376      | Protein A  |
| 5220-0377      | Goat Anti-Rat IgG (H+L)  |
| 5220-0378      | Goat Anti-Guinea Pig IgG (H+L)                                     |
| 5220-0379      | Goat Anti-Cat IgG (H+L)  |
| 5220-0380      | Goat Anti-Horse IgG (H+L)  |
| 5220-0382      | Goat Anti-Peromyscus leucopus IgG (H+L)                            |
| 5220-0383      | Goat Anti-Bovine IgG (H+L)   |
| 5220-0384      | Rabbit Anti-Goat IgG (H+L)   |
| 5220-0385      | Goat Anti-Swine IgG (H+L)  |
| 5220-0386      | Goat Anti-Mouse IgA (alpha)  |
| 5220-0387      | Goat Anti-Dog IgG (H+L)  |
| 5220-0388      | Rabbit Anti-Sheep IgG (H+L)  |
| 5220-0389      | Goat Anti-Chicken IgG (H+L)  |
| 5220-0390      | Goat Anti-Human IgG (gamma), F(ab)2 Fragment                       |
| 5220-0391      | Goat Anti-Human IgM (my), F(ab)2 Fragment                          |
| 5220-0392      | Goat Anti-Human IgG (H+L), F(ab)2 Fragment                         |
| 5220-0393      | Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed  |
| 5220-0394      | Goat Anti-Mouse IgG (gamma), F(ab)2 Fragment, Human Serum Adsorbed |
| 5220-0395      | Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed   |
| 5220-0399      | Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed   |
| 5220-0400      | Goat Anti-Human IgG (gamma), F(ab)2 Fragment                       |
| 5220-0456      | Goat Anti-Human IgG (gamma), liquid                                |
| 5220-0457      | Goat Anti-Human IgM (my), liquid                                   |
| 5220-0458      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed, liquid           |
| 5220-0459      | Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed, liquid              |
| 5220-0460      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed, liquid          |
| 5220-0462      | Goat Anti-Human IgG (gamma), liquid                                |
| 5220-0463      | Goat Anti-Human IgM (my), liquid                                   |
| 5220-0464      | Goat Anti-Human IgG (H+L), liquid                                  |
| 5220-0465      | Goat Anti-Rabbit IgG (H+L), liquid                                 |
| 5220-0466      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed, liquid           |
| 5220-0467      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed, liquid          |
| 5220-0468      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed, liquid            |
| 5230-0284      | Goat Anti-Human IgA (alpha)  |
| 5230-0285      | Goat Anti-Human IgG (gamma)  |
| 5230-0286      | Goat Anti-Human IgM (my)   |
| 5230-0287      | Goat Anti-Human IgE (epsilon)                                      |
| 5230-0288      | Goat Anti-Human IgG (H+L)  |
| 5230-0289      | Goat Anti-Human IgA+IgG+IgM (H+L)                                  |
| 5230-0290      | Goat Anti-Human IgA + IgG + IgM (H+L), Mouse Serum Adsorbed        |
| 5230-0291      | Goat Anti-Human IgG (Fc)   |
| 5230-0292      | Goat Anti-Bovine IgM (my)  |
| 5230-0293      | Goat Anti-Bovine IgG (H+L)   |
| 5230-0294      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0295      | Goat Anti-Swine IgG (gamma)  |
| 5230-0297      | Goat Anti-Swine IgG (H+L)  |
| 5230-0298      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0299      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                   |
| 5230-0301      | Goat Anti-Rat IgG (H+L)  |
| 5230-0303      | Goat Anti-Guinea Pig IgG (H+L)                                     |
| 5230-0305      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed                  |
| 5230-0307      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0308      | Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed        |
| 5230-0309      | Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed              |
| 5230-0310      | Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed              |
| 5230-0311      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5230-0312      | Goat Anti-Dog IgG (gamma)  |
| 5230-0313      | Goat Anti-Dog IgM (my)   |
| 5230-0314      | Goat Anti-Dog IgG (H+L)  |
| 5230-0315      | Goat Anti-Cat IgG (gamma)  |
| 5230-0316      | Goat Anti-Cat IgM (my)   |
| 5230-0317      | Goat Anti-Cat IgG (H+L)  |
| 5230-0318      | Goat Anti-Horse IgG (gamma)  |
| 5230-0320      | Goat Anti-Horse IgG (H+L)  |
| 5230-0322      | Rabbit Anti-Sheep IgG (gamma)                                      |
| 5230-0324      | Rabbit Anti-Sheep IgG (H+L)  |
| 5230-0325      | Goat Anti-Chicken IgG (H+L)  |
| 5230-0326      | Goat Anti-Duck IgG (H+L)   |
| 5230-0327      | Goat Anti-Turkey IgG (H+L)   |
| 5230-0329      | Rabbit Anti-Deer IgG (H+L)   |
| 5230-0332      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0336      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0338      | Goat Anti-Human IgG (H+L)  |

| Article Number | Description  |
|----------------|--|
| 5230-0339      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0340      | Goat Anti-Rat IgG (H+L)  |
| 5230-0341      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0342      | Goat Anti-Human IgG (H+L)  |
| 5230-0343      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0344      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0345      | Goat Anti-Human IgG (H+L)  |
| 5230-0347      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0348      | Goat Anti-Human IgG (H+L)  |
| 5230-0349      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0350      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0355      | Goat Anti-Human IgG (gamma)  |
| 5230-0358      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0359      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0361      | Goat Anti-Rat IgG (H+L)  |
| 5230-0362      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed                  |
| 5230-0364      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0366      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5230-0368      | Goat Anti-Human IgG (gamma)  |
| 5230-0369      | Goat Anti-Human IgM (my)   |
| 5230-0370      | Goat Anti-Human IgG (H+L)  |
| 5230-0371      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0372      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0373      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                   |
| 5230-0374      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed                  |
| 5230-0375      | Goat Anti-Mouse IgM (my), Human Serum Adsorbed                     |
| 5230-0376      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0377      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5230-0378      | Rabbit Anti-Sheep IgG (H+L)  |
| 5230-0379      | Goat Anti-Chicken IgG (H+L)  |
| 5230-0380      | Goat Anti-Human IgG (gamma)  |
| 5230-0381      | Goat Anti-Human IgM (my)   |
| 5230-0382      | Goat Anti-Human IgG (H+L)  |
| 5230-0383      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0385      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0386      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                   |
| 5230-0387      | Goat Anti-Rat IgG (H+L)  |
| 5230-0388      | Goat Anti-Guinea Pig IgG (H+L)                                     |
| 5230-0389      | Goat Anti-Mouse IgG (gamma)  |
| 5230-0390      | Goat Anti-Mouse IgM (my), Human Serum Adsorbed                     |
| 5230-0391      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0392      | Goat Anti-Mouse IgG+IgM (H+L), Human Serum Adsorbed                |
| 5230-0393      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5230-0394      | Goat Anti-Chicken IgG (H+L)  |
| 5230-0395      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0397      | Goat Anti-Guinea Pig IgG (H+L)                                     |
| 5230-0400      | Goat Anti-Human IgG (H+L)  |
| 5230-0401      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0402      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0403      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                   |
| 5230-0404      | Goat Anti-Rat IgG (H+L)  |
| 5230-0405      | Goat Anti-Guinea Pig IgG (H+L)                                     |
| 5230-0406      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0407      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5230-0409      | Goat Anti-Chicken IgG (H+L)  |
| 5230-0410      | Goat Anti-Human IgG (H+L)  |
| 5230-0411      | Rabbit Anti-Goat IgG (H+L)   |
| 5230-0412      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0413      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                   |
| 5230-0415      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0416      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed           |
| 5230-0418      | Goat Anti-Chicken IgG (H+L)  |
| 5230-0420      | Goat Anti-Mouse IgG (H+L)  |
| 5230-0421      | Goat Anti-Chicken IgG (H+L)  |
| 5230-0422      | Goat Anti-Monkey IgG (gamma)                                       |
| 5230-0423      | Goat Anti-Monkey IgM (my)  |
| 5230-0425      | Goat Anti-Human IgG (H+L)  |
| 5230-0426      | Goat Anti-Rabbit IgG (H+L)   |
| 5230-0427      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                    |
| 5230-0429      | Goat Anti-Human IgG (H+L), F(ab)2 Fragment                         |
| 5230-0431      | Goat Anti-Mouse IgG (gamma), F(ab)2 Fragment, Human Serum Adsorbed |
| 5230-0432      | Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed  |
| 5230-0436      | Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed  |
| 5230-0438      | Goat Anti-Human IgG (H+L), F(ab)2 Fragment                         |
| 5230-0440      | Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed   |
| 5230-0441      | Goat Anti-Human IgG (gamma), F(ab)2 Fragment                       |

| Article Number | Description   |
|----------------|---|
| 5230-0442      | Goat Anti-Human IgM (my), F(ab)2 Fragment                         |
| 5230-0443      | Goat Anti-Human IgG (H+L), F(ab)2 Fragment                        |
| 5230-0444      | Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed |
| 5230-0446      | Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed  |
| 5260-0027      | Goat Anti-Human IgA (alpha)                                       |
| 5260-0028      | Goat Anti-Human IgG (gamma)                                       |
| 5260-0029      | Goat Anti-Human IgM (my)  |
| 5260-0030      | Goat Anti-Human IgE (epsilon)                                     |
| 5260-0031      | Goat Anti-Human IgG (H+L)   |
| 5260-0032      | Goat Anti-Human IgA + IgG + IgM (H+L)                             |
| 5260-0033      | Goat Anti-Bovine IgM (my)   |
| 5260-0034      | Goat Anti-Bovine IgG (H+L)  |
| 5260-0035      | Rabbit Anti-Goat IgG (H+L)  |
| 5260-0036      | Goat Anti-Swine IgG (gamma)                                       |
| 5260-0037      | Goat Anti-Swine IgM (my)  |
| 5260-0038      | Goat Anti-Rabbit IgG (H+L)  |
| 5260-0039      | Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed                  |
| 5260-0041      | Goat Anti-Rat IgM (my)  |
| 5260-0042      | Goat Anti-Rat IgG (H+L)   |
| 5260-0043      | Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed                     |
| 5260-0044      | Goat Anti-Guinea Pig IgG (H+L)                                    |
| 5260-0045      | Goat Anti-Mouse IgA (alpha), Human Serum Adsorbed                 |
| 5260-0046      | Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed                 |
| 5260-0047      | Goat Anti-Mouse IgM (my), Human Serum Adsorbed                    |
| 5260-0048      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                   |
| 5260-0049      | Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed       |
| 5260-0050      | Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed             |
| 5260-0051      | Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed             |
| 5260-0052      | Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed          |
| 5260-0053      | Goat Anti-Dog IgG (H+L)   |
| 5260-0054      | Goat Anti-Cat IgG (H+L)   |
| 5260-0055      | Goat Anti-Horse IgG (gamma)                                       |
| 5260-0056      | Goat Anti-Horse IgM (my)  |
| 5260-0057      | Rabbit Anti-Sheep IgG (H+L)                                       |
| 5260-0058      | Goat Anti-Chicken IgG (H+L)                                       |
| 5260-0059      | Goat Anti-Human IgG (H+L)   |
| 5260-0060      | Goat Anti-Rabbit IgG (H+L)  |
| 5260-0061      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed                   |
| 5260-0062      | Goat Anti-Human IgG (gamma), F(ab)2 Fragment                      |
| 5260-0063      | Goat Anti-Human IgG (H+L), F(ab)2 Fragment                        |
| 5260-0064      | Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment                        |
| 5270-0017      | Streptavidin  |
| 5270-0018      | Streptavidin  |
| 5270-0019      | Streptavidin  |
| 5270-0020      | Streptavidin  |
| 5270-0021      | Streptavidin, liquid  |
| 5270-0022      | Streptavidin  |
| 5270-0023      | Streptavidin  |
| 5270-0024      | Streptavidin  |
| 5270-0025      | Streptavidin  |
| 5270-0027      | Streptavidin  |
| 5270-0028      | Streptavidin  |
| 5270-0029      | Streptavidin  |
| 5270-0030      | Streptavidin  |
| 5310-0310      | BacTrace Rabbit Anti-Shigella Species                             |
| 5310-0311      | BacTrace Rabbit Anti-Vibrio Species                               |
| 5310-0312      | BacTrace Rabbit Anti-Legionella Species                           |
| 5310-0313      | BacTrace Rabbit Anti-Yersinia Species                             |
| 5310-0314      | BacTrace Goat Anti-Staphylococcus aureus                          |
| 5310-0315      | BacTrace Goat Anti-Streptococcus pyogenes Group A                 |
| 5310-0316      | BacTrace Goat Anti-Streptococcus agalactiae Group B               |
| 5310-0317      | BacTrace Rabbit Anti-Vibrio cholerae                              |
| 5310-0318      | BacTrace Rabbit Anti-Clavibacter michiganensis                    |
| 5310-0319      | BacTrace Goat Anti-Listeria, genus-specific                       |
| 5310-0320      | BacTrace Goat Anti-Listeria, High Sensitivity                     |
| 5310-0321      | BacTrace Goat Anti-Salmonella CSA-Plus                            |
| 5310-0322      | BacTrace Goat Anti-Salmonella, CSA-1                              |
| 5310-0323      | BacTrace Goat Anti-Salmonella, CSA-1, Molecular Grade             |
| 5310-0324      | BacTrace Goat Anti-Campylobacter Species                          |
| 5310-0325      | BacTrace Goat Anti-Helicobacter pylori                            |
| 5310-0326      | BacTrace Goat Anti-E. coli O157:H7                                |
| 5310-0327      | BacTrace Goat Anti-E. coli O157:H7, Molecular Grade               |
| 5310-0328      | BacTrace Goat Anti-E. coli O111 Species                           |
| 5310-0329      | BacTrace Goat Anti-E. coli O26 Species                            |
| 5310-0330      | BacTrace Goat Anti-E. coli O103 Species                           |
| 5310-0331      | BacTrace Goat Anti-E. coli O145 Species                           |

| Article Number | Description  |
|----------------|--|
| 5310-0332      | BacTrace Goat Anti-E. coli O121 Species                        |
| 5310-0333      | BacTrace Goat Anti-E. coli O45 Species                         |
| 5310-0334      | BacTrace Goat Anti-E. coli O104 Species                        |
| 5310-0335      | BacTrace Goat Anti-E. coli O91 Species                         |
| 5310-0336      | BacTrace Goat Anti-Renibacterium salmoninarum                  |
| 5310-0337      | BacTrace Goat Anti-Borrelia burgdorferi                        |
| 5310-0338      | BacTrace Goat Anti-Borrelia Species                            |
| 5310-0339      | BacTrace Goat Anti-Staphylococcus aureus                       |
| 5310-0340      | BacTrace Goat Anti-Streptococcus pyogenes Group A              |
| 5310-0341      | BacTrace Anti-E. coli O103 Species                             |
| 5310-0342      | BacTrace Anti-E. coli O145 Species                             |
| 5310-0343      | BacTrace Anti-E. coli O121 Species                             |
| 5310-0344      | BacTrace Anti-E. coli O45 Species                              |
| 5310-0345      | BacTrace Anti-E. coli O157 Species                             |
| 5310-0346      | BacTrace Anti-E. coli O157 Species                             |
| 5310-0347      | BacTrace Anti-E. coli O111 Species                             |
| 5310-0348      | BacTrace Anti-E. coli O26 Species                              |
| 5310-0349      | BacTrace Goat Anti-Streptococcus agalactiae Group B            |
| 5310-0350      | BacTrace Rabbit Anti-Clavibacter michiganensis                 |
| 5310-0351      | BacTrace Anti-Salmonella CSA-Plus                              |
| 5310-0352      | BacTrace Goat Anti-E. coli O111 Species                        |
| 5310-0353      | BacTrace Goat Anti-E. coli O26 Species                         |
| 5310-0354      | BacTrace Goat Anti-E. coli O103 Species                        |
| 5310-0355      | BacTrace Goat Anti-E. coli O145 Species                        |
| 5310-0356      | BacTrace Goat Anti-E. coli O121 Species                        |
| 5310-0357      | BacTrace Goat Anti-E. coli O45 Species                         |
| 5310-0358      | BacTrace Goat Anti-E. coli O104 Species                        |
| 5310-0359      | BacTrace Goat Anti-E. coli O91 Species                         |
| 5320-0027      | BacTrace® Anti-Shigella Species, Peroxidase-labeled            |
| 5320-0028      | BacTrace Rabbit Anti-Vibrio Species                            |
| 5320-0029      | BacTrace® Anti-Legionella Species Antibody, Peroxidase labeled |
| 5320-0030      | BacTrace Rabbit Anti-Yersinia Species                          |
| 5320-0032      | BacTrace Goat Anti-Salmonella, CSA-1                           |
| 5320-0033      | BacTrace Goat Anti-Campylobacter Species                       |
| 5320-0034      | BacTrace Goat Anti-Helicobacter pylori                         |
| 5320-0035      | BacTrace Goat Anti-E. coli O157:H7                             |
| 5320-0036      | BacTrace Goat Anti-E. coli O111 Species                        |
| 5320-0037      | BacTrace Goat Anti-E. coli O26 Species                         |
| 5320-0038      | BacTrace Goat Anti-E. coli O103 Species                        |
| 5320-0039      | BacTrace Goat Anti-E. coli O145 Species                        |
| 5320-0040      | BacTrace Goat Anti-E. coli O121 Species                        |
| 5320-0041      | BacTrace Goat Anti-E. coli O45 Species                         |
| 5320-0042      | BacTrace Goat Anti-Renibacterium salmoninarum                  |
| 5320-0044      | BacTrace Goat Anti-Salmonella, CSA-1                           |
| 5320-0045      | BacTrace Goat Anti-Campylobacter Species                       |
| 5320-0047      | BacTrace Goat Anti-Listeria, genus-specific                    |
| 5320-0048      | BacTrace Goat Anti-Borrelia burgdorferi                        |
| 5330-0059      | BacTrace Goat Anti-Salmonella, CSA-1                           |
| 5330-0060      | BacTrace Goat Anti-Campylobacter Species                       |
| 5330-0062      | BacTrace Goat Anti-E. coli O157:H7                             |
| 5330-0063      | BacTrace Goat Anti-Renibacterium salmoninarum                  |
| 5330-0064      | BacTrace Goat Anti-Borrelia burgdorferi                        |
| 5330-0065      | BacTrace Goat Anti-Borrelia Species                            |
| 5330-0066      | BacTrace Rabbit Anti-Shigella Species                          |
| 5330-0067      | BacTrace Rabbit Anti-Vibrio Species                            |
| 5330-0068      | BacTrace Rabbit Anti-Legionella Species                        |
| 5330-0069      | BacTrace Rabbit Anti-Yersinia Species                          |
| 5330-0070      | BacTrace Goat Anti-E. coli O111 Species                        |
| 5330-0071      | BacTrace Goat Anti-E. coli O26 Species                         |
| 5330-0072      | BacTrace Goat Anti-E. coli O103 Species                        |
| 5330-0073      | BacTrace Goat Anti-E. coli O145 Species                        |
| 5330-0074      | BacTrace Goat Anti-E. coli O121 Species                        |
| 5330-0075      | BacTrace Goat Anti-E. coli O45 Species                         |
| 5350-0036      | BacTrace Goat Anti-Salmonella, CSA-1, Molecular Grade          |
| 5350-0037      | BacTrace Goat Anti-E. coli O157:H7, Molecular Grade            |
| 5350-0039      | BacTrace Goat Anti-E. coli O111 Species                        |
| 5350-0040      | BacTrace Goat Anti-E. coli O26 Species                         |
| 5350-0041      | BacTrace Goat Anti-E. coli O103 Species                        |
| 5350-0042      | BacTrace Goat Anti-E. coli O145 Species                        |
| 5350-0043      | BacTrace Goat Anti-E. coli O121 Species                        |
| 5350-0044      | BacTrace Goat Anti-E. coli O45 Species                         |
| 5360-0030      | BacTrace Goat Anti-Listeria Species                            |
| 5360-0031      | BacTrace Goat Anti-Salmonella, CSA-1                           |
| 5360-0032      | BacTrace Goat Anti-E. coli O157:H7                             |
| 5360-0033      | BacTrace Goat Anti-E. coli O111 Species                        |
| 5360-0034      | BacTrace Goat Anti-E. coli O26 Species                         |

| Article Number | Description   |
|----------------|---|
| 5360-0035      | BacTrace Goat Anti-E. coli O103 Species                 |
| 5360-0036      | BacTrace Goat Anti-E. coli O145 Species                 |
| 5360-0037      | BacTrace Goat Anti-E. coli O121 Species                 |
| 5360-0038      | BacTrace Goat Anti-E. coli O45 Species                  |
| 5450-0009      | Goat Anti-Human IgG (H+L), liquid                       |
| 5450-0010      | Goat Anti-Rabbit IgG (H+L), liquid                      |
| 5450-0011      | Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed, liquid |
| 5570-0006      | Goat Anti-Mouse IgG (H+L), Human Serum Absorbed, liquid |
| 5570-0007      | Goat Anti-Rabbit IgG (H+L), liquid                      |
| 5570-0008      | Goat Anti-Rat IgG (H+L)                                 |
| 5570-0009      | Rabbit Anti-Goat IgG (H+L)                              |
| 5950-0004      | Streptavidin, liquid                                    |
| 5950-0005      | Streptavidin, liquid                                    |

Article Number Description

---