

**Information zur Erstellung von Sicherheitsdatenblättern**

**Produkt:** alle Antikörper und Antikörper-Konjugate  
**Hersteller:** SeraCare Life Sciences  
910 Clopper Road, Gaithersburg, Maryland 20878 USA  
**Artikelnummern:** siehe Anlage

Alle in der Anlage aufgeführten Antikörper und Antikörper-Konjugate wurden gemäß der Verordnung (EG) Nr. 1907/2006 (REACH) und der Verordnung (EG) Nr. 1272/2008 (GHS/CLP) eingestuft.

Bei keinem der Antikörper und keinem der Antikörper-Konjugate handelt es sich nach Einstufung gemäß vorbenannter Verordnungen um gefährliche Produkte. Somit wurden keine Sicherheitsdatenblätter erstellt.

Beim Umgang mit den Antikörpern und Antikörper-Konjugaten sind die allgemeinen und spezifischen Sicherheitsvorkehrungen für Laboratorien zu beachten und die entsprechenden Vorschriften einzuhalten.

Ausgabedatum: Februar 2020

**medac**Produktsicherheit  
Safety, Health & Environment  
Tel.: +49 (0)4103 8006 0  
Fax: +49 (0)4103 8006 100  
E-Mail: productsafety@medac.deTheaterstrasse 6  
22880 Wedel  
GermanyGeschäftsführer:  
Jens Denker  
Jörg Hans  
Nikolaus Graf Stolberg  
Heiner WillRegistergericht:  
Pinneberg HRB 12042 PIAnlage

Artikelnummer	Beschreibung
5210-0154	Goat Anti-Human IgA (alpha)
5210-0155	Goat Anti-Human IgG (gamma)
5210-0157	Goat Anti-Human IgM (my)
5210-0158	Goat Anti-Human IgE (epsilon)
5210-0159	Goat Anti-Human IgG (H+L)
5210-0160	Goat Anti-Human IgA+IgG+IgM (H+L)
5210-0164	Goat Anti-Human IgA + IgG + IgM (H+L), Mouse Serum Adsorbed
5210-0165	Goat Anti-Human IgG (Fc)
5210-0168	Goat Anti-Bovine IgM (my)
5210-0169	Goat Anti-Bovine IgG (H+L)
5210-0170	Rabbit Anti-Goat IgG (H+L)
5210-0171	Goat Anti-Swine IgG (gamma)
5210-0172	Goat Anti-Swine IgM (my)
5210-0173	Goat Anti-Swine IgG (H+L)
5210-0174	Goat Anti-Rabbit IgG (H+L)
5210-0176	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5210-0177	Goat Anti-Rat IgG (gamma)
5210-0178	Goat Anti-Rat IgM (my)
5210-0179	Goat Anti-Rat IgG (H+L)
5210-0180	Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed
5210-0181	Goat Anti-Guinea Pig IgG (H+L)
5210-0182	Goat Anti-Mouse IgA (alpha)
5210-0183	Goat Anti-Mouse IgG (gamma)
5210-0184	Goat Anti-Mouse IgM (my), Human Serum Adsorbed
5210-0185	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5210-0187	Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed
5210-0188	Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed
5210-0191	Goat Anti-Mouse IgG (H+L), Rat Serum Adsorbed
5210-0192	Goat Anti-Mouse IgG (H+L), Rabbit Serum Adsorbed
5210-0193	Goat Anti-Dog IgG (gamma)
5210-0194	Goat Anti-Dog IgM (my)
5210-0195	Goat Anti-Dog IgG (H+L)
5210-0197	Goat Anti-Cat IgM (my)
5210-0198	Goat Anti-Cat IgG (H+L)
5210-0199	Goat Anti-Horse IgG (gamma)
5210-0200	Goat Anti-Horse IgM (my)
5210-0201	Goat Anti-Horse IgG (H+L)
5210-0203	Rabbit Anti-Sheep IgG (gamma)
5210-0204	Rabbit Anti-Sheep IgG (H+L)
5210-0205	Goat Anti-Chicken IgG (H+L)
5210-0206	Goat Anti-Duck IgG (H+L)
5210-0207	Goat Anti-Turkey IgG (H+L)
5210-0209	Rabbit Anti-Deer IgG (H+L)
5210-0210	Goat Anti-Peromyscus leucopus IgG (H+L)
5210-0211	Goat Anti-FITC
5210-0212	Goat Anti-PE
5210-0213	Goat Anti-APC
5210-0214	Goat Anti-PerCP
5210-0215	Goat Anti-Monkey IgA (alpha), liquid
5210-0216	Goat Anti-Monkey IgG (gamma), liquid
5210-0217	Goat Anti-Monkey IgM (my), liquid
5210-0297	Goat Anti-Human IgG (gamma), Fab2 Fragment
5210-0298	Goat Anti-Human IgM (my), F(ab)2 Fragment
5210-0299	Goat Anti-Human IgG (H+L), F(ab)2 Fragment
5210-0300	Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5210-0302	Goat Anti-Mouse IgG (gamma), F(ab)2 Fragment, Human Serum Adsorbed
5210-0303	Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5220-0277	Goat Anti-Human IgG (H+L)
5220-0278	Goat Anti-Human IgA+IgG+IgM (H+L), Mouse Serum Adsorbed
5220-0279	Goat Anti-Human IgG (Fc)
5220-0280	Goat Anti-Bovine IgM (my)
5220-0281	Goat Anti-Swine IgG (gamma)
5220-0282	Goat Anti-Swine IgM (my)
5220-0283	Goat Anti-Rabbit IgG (H+L)
5220-0284	Goat Anti-Rat IgG (gamma)
5220-0285	Goat Anti-Rat IgM (my)
5220-0286	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5220-0287	Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed
5220-0288	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5220-0289	Goat Anti-Dog IgG (gamma)
5220-0290	Goat Anti-Dog IgM (my)
5220-0291	Goat Anti-Cat IgG (gamma)
5220-0292	Goat Anti-Cat IgM (my)

Artikelnummer	Beschreibung
5220-0293	Goat Anti-Horse IgG (gamma)
5220-0294	Goat Anti-Horse IgM (my)
5220-0295	Rabbit Anti-Sheep IgG (gamma)
5220-0296	Goat Anti-Duck IgG (H+L)
5220-0298	Rabbit Anti-Deer IgG (H+L)
5220-0299	Goat Anti-FITC
5220-0300	Goat Anti-PE
5220-0301	Goat Anti-APC
5220-0302	Goat Anti-PerCP
5220-0303	Goat Anti-Human IgG (H+L)
5220-0304	Goat Anti-Human IgG (Fc)
5220-0305	Goat Anti-Bovine IgM (my)
5220-0306	Goat Anti-Swine IgG (gamma)
5220-0307	Goat Anti-Swine IgM (my)
5220-0308	Goat Anti-Rabbit IgG (H+L)
5220-0309	Goat Anti-Rat IgG (gamma)
5220-0310	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5220-0311	Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed
5220-0312	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5220-0314	Goat Anti-Dog IgM (my)
5220-0315	Goat Anti-Cat IgG (gamma)
5220-0316	Goat Anti-Horse IgG (gamma)
5220-0317	Goat Anti-Horse IgM (my)
5220-0318	Rabbit Anti-Sheep IgG (gamma)
5220-0320	Goat Anti-Ferret IgG (H+L)
5220-0322	Goat Anti-Dog IgG (gamma)
5220-0323	Goat Anti-FITC
5220-0324	Goat Anti-PE
5220-0325	Goat Anti-APC
5220-0326	Goat Anti-PerCP
5220-0327	Goat Anti-Human IgG (gamma)
5220-0328	Goat Anti-Human IgM (my)
5220-0329	Goat Anti-Human IgE (epsilon)
5220-0330	Goat Anti-Human IgG (H+L)
5220-0331	Goat Anti-Human IgA + IgG + IgM (H+L)
5220-0332	Goat Anti-Monkey IgA (alpha)
5220-0333	Goat Anti-Monkey IgG (gamma)
5220-0334	Goat Anti-Monkey IgM (my)
5220-0335	Goat Anti-Rabbit IgG (H+L), Extra Serum Adsorbed
5220-0336	Goat Anti-Rabbit IgG (H+L)
5220-0337	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5220-0338	Goat Anti-Mouse IgG (H+L), Extra Serum Adsorbed
5220-0339	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed
5220-0340	Goat Anti-Mouse IgM (my), Human Serum Adsorbed
5220-0341	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5220-0342	Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed
5220-0343	Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed
5220-0347	Goat Anti-Human IgA (alpha)
5220-0348	Goat Anti-Human IgG (gamma)
5220-0349	Goat Anti-Human IgM (my)
5220-0350	Goat Anti-Human IgE (epsilon)
5220-0351	Goat Anti-Human IgG (H+L)
5220-0352	Goat Anti-Human IgA + IgG + IgM (H+L)
5220-0353	Goat Anti-Rabbit IgG (H+L)
5220-0354	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5220-0355	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed
5220-0356	Goat Anti-Mouse IgM (my), Human Serum Adsorbed
5220-0357	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5220-0358	Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed
5220-0359	Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed
5220-0360	Goat Anti-Human IgA (alpha)
5220-0361	Goat Anti-Bovine IgG (H+L)
5220-0362	Rabbit Anti-Goat IgG (H+L)
5220-0363	Goat Anti-Swine IgG (H+L)
5220-0364	Goat Anti-Rat IgG (H+L)
5220-0365	Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed
5220-0366	Goat Anti-Guinea Pig IgG (H+L)
5220-0367	Goat Anti-Mouse IgA (alpha), Human Serum Adsorbed
5220-0368	Goat Anti-Dog IgG (H+L)
5220-0369	Goat Anti-Cat IgG (H+L)
5220-0370	Goat Anti-Horse IgG (H+L)
5220-0371	Goat Anti-Hamster IgG (H+L)
5220-0372	Rabbit Anti-Sheep IgG (H+L)

Artikelnummer	Beschreibung
5220-0373	Goat Anti-Chicken IgG (H+L)
5220-0374	Goat Anti-Turkey IgG (H+L)
5220-0375	Goat Anti-Peromyscus leucopus IgG (H+L)
5220-0376	Protein A
5220-0377	Goat Anti-Rat IgG (H+L)
5220-0378	Goat Anti-Guinea Pig IgG (H+L)
5220-0379	Goat Anti-Cat IgG (H+L)
5220-0380	Goat Anti-Horse IgG (H+L)
5220-0382	Goat Anti-Peromyscus leucopus IgG (H+L)
5220-0383	Goat Anti-Bovine IgG (H+L)
5220-0384	Rabbit Anti-Goat IgG (H+L)
5220-0385	Goat Anti-Swine IgG (H+L)
5220-0386	Goat Anti-Mouse IgA (alpha)
5220-0387	Goat Anti-Dog IgG (H+L)
5220-0388	Rabbit Anti-Sheep IgG (H+L)
5220-0389	Goat Anti-Chicken IgG (H+L)
5220-0390	Goat Anti-Human IgG (gamma), F(ab)2 Fragment
5220-0391	Goat Anti-Human IgM (my), F(ab)2 Fragment
5220-0392	Goat Anti-Human IgG (H+L), F(ab)2 Fragment
5220-0393	Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5220-0394	Goat Anti-Mouse IgG (gamma), F(ab)2 Fragment, Human Serum Adsorbed
5220-0395	Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5220-0399	Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5220-0400	Goat Anti-Human IgG (gamma), F(ab)2 Fragment
5220-0456	Goat Anti-Human IgG (gamma), liquid
5220-0457	Goat Anti-Human IgM (my), liquid
5220-0458	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed, liquid
5220-0459	Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed, liquid
5220-0460	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed, liquid
5220-0462	Goat Anti-Human IgG (gamma), liquid
5220-0463	Goat Anti-Human IgM (my), liquid
5220-0464	Goat Anti-Human IgG (H+L), liquid
5220-0465	Goat Anti-Rabbit IgG (H+L), liquid
5220-0466	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed, liquid
5220-0467	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed, liquid
5220-0468	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed, liquid
5230-0284	Goat Anti-Human IgA (alpha)
5230-0285	Goat Anti-Human IgG (gamma)
5230-0286	Goat Anti-Human IgM (my)
5230-0287	Goat Anti-Human IgE (epsilon)
5230-0288	Goat Anti-Human IgG (H+L)
5230-0289	Goat Anti-Human IgA+IgG+IgM (H+L)
5230-0290	Goat Anti-Human IgA + IgG + IgM (H+L), Mouse Serum Adsorbed
5230-0291	Goat Anti-Human IgG (Fc)
5230-0292	Goat Anti-Bovine IgM (my)
5230-0293	Goat Anti-Bovine IgG (H+L)
5230-0294	Rabbit Anti-Goat IgG (H+L)
5230-0295	Goat Anti-Swine IgG (gamma)
5230-0297	Goat Anti-Swine IgG (H+L)
5230-0298	Goat Anti-Rabbit IgG (H+L)
5230-0299	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5230-0301	Goat Anti-Rat IgG (H+L)
5230-0303	Goat Anti-Guinea Pig IgG (H+L)
5230-0305	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed
5230-0307	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0308	Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed
5230-0309	Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed
5230-0310	Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed
5230-0311	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5230-0312	Goat Anti-Dog IgG (gamma)
5230-0313	Goat Anti-Dog IgM (my)
5230-0314	Goat Anti-Dog IgG (H+L)
5230-0315	Goat Anti-Cat IgG (gamma)
5230-0316	Goat Anti-Cat IgM (my)
5230-0317	Goat Anti-Cat IgG (H+L)
5230-0318	Goat Anti-Horse IgG (gamma)
5230-0320	Goat Anti-Horse IgG (H+L)
5230-0322	Rabbit Anti-Sheep IgG (gamma)
5230-0324	Rabbit Anti-Sheep IgG (H+L)
5230-0325	Goat Anti-Chicken IgG (H+L)
5230-0326	Goat Anti-Duck IgG (H+L)
5230-0327	Goat Anti-Turkey IgG (H+L)
5230-0329	Rabbit Anti-Deer IgG (H+L)
5230-0332	Goat Anti-Rabbit IgG (H+L)
5230-0336	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0338	Goat Anti-Human IgG (H+L)

Artikelnummer	Beschreibung
5230-0339	Goat Anti-Rabbit IgG (H+L)
5230-0340	Goat Anti-Rat IgG (H+L)
5230-0341	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0342	Goat Anti-Human IgG (H+L)
5230-0343	Goat Anti-Rabbit IgG (H+L)
5230-0344	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0345	Goat Anti-Human IgG (H+L)
5230-0347	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0348	Goat Anti-Human IgG (H+L)
5230-0349	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0350	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0355	Goat Anti-Human IgG (gamma)
5230-0358	Rabbit Anti-Goat IgG (H+L)
5230-0359	Goat Anti-Rabbit IgG (H+L)
5230-0361	Goat Anti-Rat IgG (H+L)
5230-0362	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed
5230-0364	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0366	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5230-0368	Goat Anti-Human IgG (gamma)
5230-0369	Goat Anti-Human IgM (my)
5230-0370	Goat Anti-Human IgG (H+L)
5230-0371	Rabbit Anti-Goat IgG (H+L)
5230-0372	Goat Anti-Rabbit IgG (H+L)
5230-0373	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5230-0374	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed
5230-0375	Goat Anti-Mouse IgM (my), Human Serum Adsorbed
5230-0376	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0377	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5230-0378	Rabbit Anti-Sheep IgG (H+L)
5230-0379	Goat Anti-Chicken IgG (H+L)
5230-0380	Goat Anti-Human IgG (gamma)
5230-0381	Goat Anti-Human IgM (my)
5230-0382	Goat Anti-Human IgG (H+L)
5230-0383	Rabbit Anti-Goat IgG (H+L)
5230-0385	Goat Anti-Rabbit IgG (H+L)
5230-0386	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5230-0387	Goat Anti-Rat IgG (H+L)
5230-0388	Goat Anti-Guinea Pig IgG (H+L)
5230-0389	Goat Anti-Mouse IgG (gamma)
5230-0390	Goat Anti-Mouse IgM (my), Human Serum Adsorbed
5230-0391	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0392	Goat Anti-Mouse IgG+IgM (H+L), Human Serum Adsorbed
5230-0393	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5230-0394	Goat Anti-Chicken IgG (H+L)
5230-0395	Rabbit Anti-Goat IgG (H+L)
5230-0397	Goat Anti-Guinea Pig IgG (H+L)
5230-0400	Goat Anti-Human IgG (H+L)
5230-0401	Rabbit Anti-Goat IgG (H+L)
5230-0402	Goat Anti-Rabbit IgG (H+L)
5230-0403	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5230-0404	Goat Anti-Rat IgG (H+L)
5230-0405	Goat Anti-Guinea Pig IgG (H+L)
5230-0406	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0407	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5230-0409	Goat Anti-Chicken IgG (H+L)
5230-0410	Goat Anti-Human IgG (H+L)
5230-0411	Rabbit Anti-Goat IgG (H+L)
5230-0412	Goat Anti-Rabbit IgG (H+L)
5230-0413	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5230-0415	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0416	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5230-0418	Goat Anti-Chicken IgG (H+L)
5230-0420	Goat Anti-Mouse IgG (H+L)
5230-0421	Goat Anti-Chicken IgG (H+L)
5230-0422	Goat Anti-Monkey IgG (gamma)
5230-0423	Goat Anti-Monkey IgM (my)
5230-0425	Goat Anti-Human IgG (H+L)
5230-0426	Goat Anti-Rabbit IgG (H+L)
5230-0427	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5230-0429	Goat Anti-Human IgG (H+L), F(ab)2 Fragment
5230-0431	Goat Anti-Mouse IgG (gamma), F(ab)2 Fragment, Human Serum Adsorbed
5230-0432	Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5230-0436	Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5230-0438	Goat Anti-Human IgG (H+L), F(ab)2 Fragment
5230-0440	Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5230-0441	Goat Anti-Human IgG (gamma), F(ab)2 Fragment

Artikelnummer	Beschreibung
5230-0442	Goat Anti-Human IgM (my), F(ab)2 Fragment
5230-0443	Goat Anti-Human IgG (H+L), F(ab)2 Fragment
5230-0444	Goat Anti-Rabbit IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5230-0446	Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment, Human Serum Adsorbed
5260-0027	Goat Anti-Human IgA (alpha)
5260-0028	Goat Anti-Human IgG (gamma)
5260-0029	Goat Anti-Human IgM (my)
5260-0030	Goat Anti-Human IgE (epsilon)
5260-0031	Goat Anti-Human IgG (H+L)
5260-0032	Goat Anti-Human IgA + IgG + IgM (H+L)
5260-0033	Goat Anti-Bovine IgM (my)
5260-0034	Goat Anti-Bovine IgG (H+L)
5260-0035	Rabbit Anti-Goat IgG (H+L)
5260-0036	Goat Anti-Swine IgG (gamma)
5260-0037	Goat Anti-Swine IgM (my)
5260-0038	Goat Anti-Rabbit IgG (H+L)
5260-0039	Goat Anti-Rabbit IgG (H+L), Human Serum Adsorbed
5260-0041	Goat Anti-Rat IgM (my)
5260-0042	Goat Anti-Rat IgG (H+L)
5260-0043	Goat Anti-Rat IgG (H+L), Mouse Serum Adsorbed
5260-0044	Goat Anti-Guinea Pig IgG (H+L)
5260-0045	Goat Anti-Mouse IgA (alpha), Human Serum Adsorbed
5260-0046	Goat Anti-Mouse IgG (gamma), Human Serum Adsorbed
5260-0047	Goat Anti-Mouse IgM (my), Human Serum Adsorbed
5260-0048	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5260-0049	Goat Anti-Mouse IgA + IgG + IgM (H+L), Human Serum Adsorbed
5260-0050	Goat Anti-Mouse IgG + IgM (H+L), Human Serum Adsorbed
5260-0051	Goat Anti-Mouse IgG (H+L), Rat + Human Serum Adsorbed
5260-0052	Goat Anti-Mouse IgG (H+L), Rabbit + Human Serum Adsorbed
5260-0053	Goat Anti-Dog IgG (H+L)
5260-0054	Goat Anti-Cat IgG (H+L)
5260-0055	Goat Anti-Horse IgG (gamma)
5260-0056	Goat Anti-Horse IgM (my)
5260-0057	Rabbit Anti-Sheep IgG (H+L)
5260-0058	Goat Anti-Chicken IgG (H+L)
5260-0059	Goat Anti-Human IgG (H+L)
5260-0060	Goat Anti-Rabbit IgG (H+L)
5260-0061	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed
5260-0062	Goat Anti-Human IgG (gamma), F(ab)2 Fragment
5260-0063	Goat Anti-Human IgG (H+L), F(ab)2 Fragment
5260-0064	Goat Anti-Mouse IgG (H+L), F(ab)2 Fragment
5270-0017	Streptavidin
5270-0018	Streptavidin
5270-0019	Streptavidin
5270-0020	Streptavidin
5270-0021	Streptavidin, liquid
5270-0022	Streptavidin
5270-0023	Streptavidin
5270-0024	Streptavidin
5270-0025	Streptavidin
5270-0027	Streptavidin
5270-0028	Streptavidin
5270-0029	Streptavidin
5270-0030	Streptavidin
5310-0310	BacTrace Rabbit Anti-Shigella Species
5310-0311	BacTrace Rabbit Anti-Vibrio Species
5310-0312	BacTrace Rabbit Anti-Legionella Species
5310-0313	BacTrace Rabbit Anti-Yersinia Species
5310-0314	BacTrace Goat Anti-Staphylococcus aureus
5310-0315	BacTrace Goat Anti-Streptococcus pyogenes Group A
5310-0316	BacTrace Goat Anti-Streptococcus agalactiae Group B
5310-0317	BacTrace Rabbit Anti-Vibrio cholerae
5310-0318	BacTrace Rabbit Anti-Clavibacter michiganensis
5310-0319	BacTrace Goat Anti-Listeria, genus-specific
5310-0320	BacTrace Goat Anti-Listeria, High Sensitivity
5310-0321	BacTrace Goat Anti-Salmonella CSA-Plus
5310-0322	BacTrace Goat Anti-Salmonella, CSA-1
5310-0323	BacTrace Goat Anti-Salmonella, CSA-1, Molecular Grade
5310-0324	BacTrace Goat Anti-Campylobacter Species
5310-0325	BacTrace Goat Anti-Helicobacter pylori
5310-0326	BacTrace Goat Anti-E. coli O157:H7
5310-0327	BacTrace Goat Anti-E. coli O157:H7, Molecular Grade
5310-0328	BacTrace Goat Anti-E. coli O111 Species
5310-0329	BacTrace Goat Anti-E. coli O26 Species
5310-0330	BacTrace Goat Anti-E. coli O103 Species
5310-0331	BacTrace Goat Anti-E. coli O145 Species

Artikelnummer	Beschreibung
5310-0332	BacTrace Goat Anti-E. coli O121 Species
5310-0333	BacTrace Goat Anti-E. coli O45 Species
5310-0334	BacTrace Goat Anti-E. coli O104 Species
5310-0335	BacTrace Goat Anti-E. coli O91 Species
5310-0336	BacTrace Goat Anti-Renibacterium salmoninarum
5310-0337	BacTrace Goat Anti-Borrelia burgdorferi
5310-0338	BacTrace Goat Anti-Borrelia Species
5310-0339	BacTrace Goat Anti-Staphylococcus aureus
5310-0340	BacTrace Goat Anti-Streptococcus pyogenes Group A
5310-0341	BacTrace Anti-E. coli O103 Species
5310-0342	BacTrace Anti-E. coli O145 Species
5310-0343	BacTrace Anti-E. coli O121 Species
5310-0344	BacTrace Anti-E. coli O45 Species
5310-0345	BacTrace Goat Anti-Salmonella CSA-1
5310-0346	BacTrace Anti-E. coli O157 Species
5310-0347	BacTrace Anti-E. coli O111 Species
5310-0348	BacTrace Anti-E. coli O26 Species
5310-0349	BacTrace Goat Anti-Streptococcus agalactiae Group B
5310-0350	BacTrace Rabbit Anti-Clavibacter michiganensis
5310-0351	BacTrace Anti-Salmonella CSA-Plus
5310-0352	BacTrace Goat Anti-E. coli O111 Species
5310-0353	BacTrace Goat Anti-E. coli O26 Species
5310-0354	BacTrace Goat Anti-E. coli O103 Species
5310-0355	BacTrace Goat Anti-E. coli O145 Species
5310-0356	BacTrace Goat Anti-E. coli O121 Species
5310-0357	BacTrace Goat Anti-E. coli O45 Species
5310-0358	BacTrace Goat Anti-E. coli O104 Species
5310-0359	BacTrace Goat Anti-E. coli O91 Species
5320-0027	BacTrace® Anti-Shigella Species, Peroxidase-labeled
5320-0028	BacTrace Rabbit Anti-Vibrio Species
5320-0029	BacTrace® Anti-Legionella Species Antibody, Peroxidase labeled
5320-0030	BacTrace Rabbit Anti-Yersinia Species
5320-0032	BacTrace Goat Anti-Salmonella, CSA-1
5320-0033	BacTrace Goat Anti-Campylobacter Species
5320-0034	BacTrace Goat Anti-Helicobacter pylori
5320-0035	BacTrace Goat Anti-E. coli O157:H7
5320-0036	BacTrace Goat Anti-E. coli O111 Species
5320-0037	BacTrace Goat Anti-E. coli O26 Species
5320-0038	BacTrace Goat Anti-E. coli O103 Species
5320-0039	BacTrace Goat Anti-E. coli O145 Species
5320-0040	BacTrace Goat Anti-E. coli O121 Species
5320-0041	BacTrace Goat Anti-E. coli O45 Species
5320-0042	BacTrace Goat Anti-Renibacterium salmoninarum
5320-0044	BacTrace Goat Anti-Salmonella, CSA-1
5320-0045	BacTrace Goat Anti-Campylobacter Species
5320-0047	BacTrace Goat Anti-Listeria, genus-specific
5320-0048	BacTrace Goat Anti-Borrelia burgdorferi
5330-0059	BacTrace Goat Anti-Salmonella, CSA-1
5330-0060	BacTrace Goat Anti-Campylobacter Species
5330-0062	BacTrace Goat Anti-E. coli O157:H7
5330-0063	BacTrace Goat Anti-Renibacterium salmoninarum
5330-0064	BacTrace Goat Anti-Borrelia burgdorferi
5330-0065	BacTrace Goat Anti-Borrelia Species
5330-0066	BacTrace Rabbit Anti-Shigella Species
5330-0067	BacTrace Rabbit Anti-Vibrio Species
5330-0068	BacTrace Rabbit Anti-Legionella Species
5330-0069	BacTrace Rabbit Anti-Yersinia Species
5330-0070	BacTrace Goat Anti-E. coli O111 Species
5330-0071	BacTrace Goat Anti-E. coli O26 Species
5330-0072	BacTrace Goat Anti-E. coli O103 Species
5330-0073	BacTrace Goat Anti-E. coli O145 Species
5330-0074	BacTrace Goat Anti-E. coli O121 Species
5330-0075	BacTrace Goat Anti-E. coli O45 Species
5350-0036	BacTrace Goat Anti-Salmonella, CSA-1, Molecular Grade
5350-0037	BacTrace Goat Anti-E. coli O157:H7, Molecular Grade
5350-0039	BacTrace Goat Anti-E. coli O111 Species
5350-0040	BacTrace Goat Anti-E. coli O26 Species
5350-0041	BacTrace Goat Anti-E. coli O103 Species
5350-0042	BacTrace Goat Anti-E. coli O145 Species
5350-0043	BacTrace Goat Anti-E. coli O121 Species
5350-0044	BacTrace Goat Anti-E. coli O45 Species
5360-0030	BacTrace Goat Anti-Listeria Species
5360-0031	BacTrace Goat Anti-Salmonella, CSA-1
5360-0032	BacTrace Goat Anti-E. coli O157:H7
5360-0033	BacTrace Goat Anti-E. coli O111 Species
5360-0034	BacTrace Goat Anti-E. coli O26 Species

Artikelnummer	Beschreibung
5360-0035	BacTrace Goat Anti-E. coli O103 Species
5360-0036	BacTrace Goat Anti-E. coli O145 Species
5360-0037	BacTrace Goat Anti-E. coli O121 Species
5360-0038	BacTrace Goat Anti-E. coli O45 Species
5450-0009	Goat Anti-Human IgG (H+L), liquid
5450-0010	Goat Anti-Rabbit IgG (H+L), liquid
5450-0011	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed, liquid
5570-0006	Goat Anti-Mouse IgG (H+L), Human Serum Adsorbed, liquid
5570-0007	Goat Anti-Rabbit IgG (H+L), liquid
5570-0008	Goat Anti-Rat IgG (H+L)
5570-0009	Rabbit Anti-Goat IgG (H+L)
5950-0004	Streptavidin, liquid
5950-0005	Streptavidin, liquid

Artikelnummer Beschreibung

---